

# Medicut Scissors

## Medicon Instruments: Made in Germany High quality, durable, long lasting

*made in Germany*

All Medicon instruments are 'Made in Germany' - according to strict standards to ensure the instruments produced are high quality and durable for long lasting use by surgeons.

Medicon has been designing and manufacturing instruments since 1941 and is now a co-operative comprising 16 member companies, of which each specialises in a product range.

Medicon instruments:

- Have both the mechanical and product engineering experience to develop modern instruments to meet surgeon's needs
- High, permanent quality & consistent designs
- Produced in Medicon facilities in Germany
- Focused on customer-oriented solutions
- Have a lifetime replacement warranty on materials and workmanship

## Medicon Medicut Scissors: ultimate cutting efficiency

- Unique ground edges ensure ultimate cutting efficiency
- Special micro serration and ground edges ensure precise, fine and minimally traumatising cuts
- Available in a wide range of sizes, lengths and designs
- Range includes scissors for delicate and extremely delicate dissecting, rhinoplasty & micro scissors
- Solid tungsten carbide (TC) alloy on the Magic cut, Supercut and TC scissors guarantee excellent cutting quality and long life-span

Medicut Scissors are perfect for plastic surgeon's who want a precise, fine cut and know that their scissors will remain sharp and accurate even after long lasting use



**medicon®**

### CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

[info@nationalsurgical.com.au](mailto:info@nationalsurgical.com.au)

[www.nationalsurgical.com.au](http://www.nationalsurgical.com.au)

Contact your National Surgical Corporation  
Account Manager for more information



BE SUPPORTED BY AN  
AUSTRALIAN COMPANY  
PROVIDING VALUE BEYOND THE PRODUCT

**NATIONAL SURGICAL CORPORATION**

*A proudly founded and owned Australian company - providing value to our customers beyond the product*

# Medicut Scissors

## magic cut

Black finish, Golden screw

- Knife edges & micro serration for a fine, minimally traumatising cut
- Welded tungsten carbide cutting edges: seamless joints, corrosion free
- Reduction of transversal cut: better visibility & precise dissection
- Ceramic coating: long lasting cutting ability



## supercut plus

Golden screw, golden rings

- Knife edges & micro serration for a fine, minimally traumatising cut
- Welded tungsten carbide cutting edges: seamless joints, corrosion free
- Reduction of transversal cut: better visibility & precise dissection



## tungsten carbide "TC"

Golden rings

- Welded tungsten carbide (TC) cutting edges: seamless joints, corrosion free



## black line

Golden screw, black rings

- Knife edges & micro serration for a fine, minimally traumatising cut
- Reduction of transversal cut: better visibility & precise dissection



## supercut

Golden screw

- Knife edges & micro serration for a fine, minimally traumatising cut



Each of the above styles are available in a large range of scissors

Check out our extensive range of scissors available for a variety of plastic surgery requirements



National Surgical Corporation  
is certified to  
AS/NZS ISO 9001:2008 (NSW & WA)



RE-SUPPORTED BY AN  
AUSTRALIAN COMPANY  
PROVIDING VALUE BEYOND THE PRODUCT

**NATIONAL SURGICAL CORPORATION**

*A proudly founded and owned Australian company - providing value to our customers beyond the product*