



# Privac-High Wound Drainage System

## *National Surgical Corporation Privac-High Wound Drainage Systems: cost effective, high suction, closed system*

### **HIGH SUCTION**

Maximum suction to 700mmHg/ 93 kPa

### **VERSATILE DRAIN FOR ANY SURGERY**

For use in any surgery where a high vacuum drain is preferred. Up to 600ml capacity

### **CLOSED SYSTEM**

Complies with Australian Infection Control Standards and OH&S guidelines

### **COMPLETE KITS AVAILABLE**

Comes with universal connection, 10FG, 14FG or 18FG drain tube, trocar, leur lock and reservoir bottle (200ml, 400ml or 600ml)

### **STERILE REPLACEMENT BOTTLES**

Sterile, single use 200ml, 400ml or 600ml reservoir bottles can be attached to the system as existing bottles fill

### **PATIENT COMFORT & MOBILITY**

Hanger allows attachment of system close to patient, on clothing or mobility aid

### **OTHER WOUND DRAINAGE OPTIONS**

Privac-Mini, Privac-Low, Privac-Tvac and Privac-Reinfusion wound drainage systems and Surgiblake silicone channel drains are also available



National Surgical Corporation  
is certified to  
AS/NZS ISO 9001:2015



# Privac-High Wound Drainage System



Product code	Description	Reservoir size	Connection	Drain size	Unit of Purchase
P24741	Privac-High set - drain, trocar, reservoir, leur lock	200ml	Universal	10 FG	Box 60
P24621		400ml	Universal	10 FG	Box 50
P21949	Privac-High OT system - reservoir, connecting tube, leur lock	400ml	Universal	-	Box 50
P24760	Privac-High set - drain, trocar, reservoir, leur lock	600ml	Universal	10 FG	Box 35
P24762		600ml	Universal	14 FG	Box 35
P24764		600ml	Universal	18 FG	Box 35
P24739	Privac-High reservoir only with leur lock	200ml	Universal	-	Box 60
P22406		600ml		-	Box 35
P21948	Privac-High ward system - reservoir only with leur lock	400ml		-	Box 50
P22407	Privac-High reservoir with connecting tube, leur lock & universal connector	600ml	Universal	-	Box 35

**National Surgical Corporation has a full range of wound drainage systems including Privac-Mini, Privac-Low, Privac-Tvac, Privac-Reinfusion and Surgiblake silicone channel drains - please refer to their brochures for more information**

## Education

National Surgical Corporation (NSC) provides trained educators to inservice all necessary staff within the hospital and home care unit in the correct operation of NSC wound drainage systems. To find out about the range of education resources available, please contact us.

### CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au

**Contact your National Surgical Corporation  
Account Manager for more information**



**NATIONAL SURGICAL CORPORATION**

**A proudly founded and owned Australian company - providing value to our customers beyond the product**