


**NEW & IMPROVED
ONE PIECE EXTRUSION**

**SOFT AND FLEXIBLE
FOR PATIENT
COMFORT**



Surgiblake Wound Drains

***National Surgical Corporation Surgiblake Wound Drains:
new and improved design, cost effective, flexible, high quality***

SUPERIOR DRAINAGE

Four flow channels for superior tissue contact, no holes reduces tissue ingrowth and less risk of blockage

SILICONE TUBING - PATIENT COMFORT

Flexible and soft silicone ensures less tissue trauma on removal of drain, making it more comfortable for the patient

SINGLE PIECE EXTRUSION

One piece extrusion reduces risk of breakage compared with two piece silicone drain designs

BENDABLE TROCAR

Trocar with 2 bendable sections to easily manipulate to desired angle

RADIO OPAQUE

To view using an x-ray if required

VARIOUS KITS AVAILABLE

Available as a round or flat drain, with or without trocar and connector to suit most wound drain systems

USE WITH ANY LOW VACUUM SYSTEM

Can be used with any low vacuum system including Privac-Low or Privac-Tvac wound drainage systems

AVAILABLE WITH 100ML BULB

Available with an easy squeeze 100ml bulb with non return valve for minor surgery

WOUND DRAINAGE OPTIONS


Privac-Mini, Privac-Low, Privac-High, Privac-Reinfusion and Privac-Tvac wound drainage systems are also available




National Surgical Corporation
is certified to
AS/NZS ISO 9001:2015



Surgiblake Wound Drains



Surgiblake
bendable trocar



Round
Surgiblake drain



Flat Surgiblake
drain



100ml Bulb as
included in SB100-02

Product Code	Description	Size	Drain size	Drain	UOP
SB500-51	Surgiblake silicone wound drain with trocar	3 x 390 x 800mm	10 FG		Box 10
SB500-52		5 x 390 x 800mm	15 FG		Box 10
SB500-53		7 x 390 x 800mm	19 FG		Box 10
SB500-53ST		7 x 390 x 800mm small trocar	19 FG		Box 10
SB500-54		8 x 390 x 800mm	Flat		Box 10
SB500-55		11 x 390 x 800mm	Flat		Box 10
SB500-41	Surgiblake silicone wound drain without trocar	3 x 390 x 800mm	10 FG		Box 10
SB500-42		5 x 390 x 800mm	15 FG		Box 10
SB500-43		11 x 390 x 800mm	19 FG		Box 10
SB500-46		8 x 350 x 800mm	24 FG		Box 10
SB500-44		8 x 390 x 800mm	Flat		Box 10
SB500-45		11 x 390 x 800mm	Flat		Box 10
SB100-02	Surgiblake silicone drain with 100ml bulb, drain and trocar	3 x 380mm	10 FG		Box 10
SB500-37	Surgiblake hubless silicone drain without trocar	7mm	Flat		Box 10
SB500-38		11mm	Flat		Box 10
SB100-12	Silicone bulb only - 100ml with non return valve	-	-		Box 10
SB200-10	Surgiblake bendable trocar (no drain)		10 FG		Box 10
SB200-15			15 FG		Box 10

CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au

Contact your National Surgical Corporation
Account Manager for more information



NATIONAL SURGICAL CORPORATION

A proudly founded and owned Australian company - providing value to our customers beyond the product