

Hermes

Closed Aluminium
Transport Case Carts



NATIONAL SURGICAL CORPORATION

Australian distributor of the Pegasus high density storage range including Hermes case carts



**TRANSPORTING EQUIPMENT INSIDE
& BETWEEN HOSPITAL DEPARTMENTS**

SECURE TRANSPORTATION

Rigid heavy-duty aluminium structure carts for transporting of sterile and non-sterile medical equipment between campuses or between central sterilisation and the receiving departments




Upper all-around bumpers for transporting inside trucks



All-around lower bumpers



Dedicated hooking structure for safe towing



COMBINATION OF LIGHT-WEIGHT AND RIGID STRUCTURE WITH CLEANING & WASHING RESISTANCE

ERGONOMIC HANDLES

Ergonomic push handles, vertical or horizontal

DOOR LOCKS

Doors can be locked with key or a seal locks

270° OPENING

Doors can be opened to 270° and locked with a magnet to the cart's outside walls

HEAVY-DUTY WHEELS

150mm heavy duty castors with or without locking

BUMPERS

All around rubber upper & bottom bumpers

TILTED TRAYS


Baskets and trays can be tilted while pulled out

PANELS

Side panels are made of durable nylon, washable and removable (SEE NEXT PAGE)

TOWING

Optional traction connection for safe towing



PANELS CAN BE REMOVED EASILY FOR WASHING

Side and central panels can be removed easily for washing using fast connectors (no tools needed)



Release locking knobs (Tool free) → Pull out middle wall (Light weight wall) → Release knobs and pull out side walls → Ready for washing (Machine wash)

NYLON SIDE PANELS

Nylon side panels are resistant to heat up to 110° degrees



CLEANING INSTRUCTIONS



- Carts are constructed of anodized light alloy and suitable for washing and chemical disinfection.
- Thermal disinfection is not recommended!
- Considering the significant number of different types of cleaning products in the market, we suggest consulting with Pegasus first.



SINGLE U CART



DOUBLE UU CART



Let NSC support you



What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.

Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has an extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!

Our range includes



Wound drains



Electrosurgical range



OT and general consumables



Surgical instruments



Laparoscopic instruments



Medichill ice packs



Protection products



Cleaning & sterilisation



Shelving & storage



High density storage

NSC has a range of shelving and storage solutions including chrome wire shelving, trolleys, compactors, high density storage, carts and height adjustable tables.

Contact us to arrange a site visit today to discuss your needs

CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au

Contact your National Surgical Corporation Account Manager for more information



BE SUPPORTED BY AN
AUSTRALIAN COMPANY
PROVIDING VALUE BEYOND THE PRODUCT



NATIONAL SURGICAL CORPORATION

A proudly founded and owned Australian company - providing value to our customers beyond the product