

THE SALUS OPEN TRANSPORTATION CART IS DESIGNED TO PROVIDE A LIGHT WEIGHT DURABLE CART FOR TRANSPORTATION OF MEDICAL SUPPLIES, CONSUMABLES, STERILE TRAYS BETWEEN DEPARTMENTS

- ✓ **LIGHTWEIGHT**
- ✓ **CLEANABLE**
- ✓ **EASY TO ASSEMBLE**
- ✓ **NO RUST**
- ✓ **NO FUTURE PAINT PEELING**

- Frame made of high-quality anodised aluminium
- Rails made from combined nylon - aluminium
- Versatile: use with ABS baskets (10cm, 20cm), ABS trays (5cm, 10cm, 20cm), wire instrument trays, flat shelves
- Easily clean and disinfect
- Various configurations: single or double cart in U or E styles in different heights

Hospital theatre • CSSD • Set up trolley • Store rooms • Emergency Dept • Wards • Day Surgeries • Clinics



1130mm height

1640mm height

1880mm height

1370mm height double

Height	Product Code U single 570 x 690mm (WxD)	Product Code E single 770 x 490mm (WxD)	Product Code U double 1050 x 690mm (WxD)	Product Code E double 1450 x 490mm (WxD)
1130mm	SLS-U-1100	SLS-E-1100	SLS-UU-1100	SLS-EE-1100
1370mm	SLS-U-1340	SLS-E-1340	SLS-UU-1340	SLS-EE-1340
1640mm	SLS-U-1580	SLS-E-1580	SLS-UU-1580	SLS-EE-1580
1880mm	SLS-U-1820	SLS-E-1820	SLS-UU-1820	SLS-EE-1820

CUSTOMER SERVICE

Phone: 1800 138 138 (free call)
info@nationalsurgical.com.au
www.nationalsurgical.com.au

Contact your National Surgical Corporation
Account Manager for more information



NATIONAL SURGICAL CORPORATION

A proudly founded and owned Australian company - providing value to our customers beyond the product