


ONE-PIECE EXTRUSION

SOFT AND FLEXIBLE  
FOR PATIENT  
COMFORT



# Surgiblake Wound Drains

*National Surgical Corporation Surgiblake Wound Drains:  
single piece extrusion, cost effective, flexible, high quality*

## **SUPERIOR DRAINAGE**

Four flow channels for superior tissue contact, no holes reduces tissue ingrowth and less risk of blockage

## **SILICONE TUBING - PATIENT COMFORT**

Flexible and soft silicone ensures less tissue trauma on removal of drain, making it more comfortable for the patient

## **SINGLE PIECE EXTRUSION**

One piece extrusion reduces risk of breakage compared with two piece silicone drain designs

## **BENDABLE TROCAR**

Trocar with 2 bendable sections to easily manipulate to desired angle

## **RADIO OPAQUE**

To view using an x-ray if required

## **VARIOUS KITS AVAILABLE**

Available as a round or flat drain, with or without trocar and connector to suit most wound drain systems

## **USE WITH ANY LOW VACUUM SYSTEM**

Can be used with any low vacuum system including Privac-Low or Privac-Tvac wound drainage systems

## **AVAILABLE WITH 100ML/400ML BULB**

Available with an easy squeeze 100ml or 400ml bulb with non-return valve for minor surgery

## **WOUND DRAINAGE OPTIONS**



























Privac-Mini, Privac-Low, Privac-High, Privac-Reinfusion and Privac-Tvac wound drainage systems are also available



National Surgical Corporation  
is certified to  
AS/NZS ISO 9001:2015



# Surgiblake Wound Drains

Product Code	Size	Drain size	Drain type	Image	UOP
Surgiblake silicone wound drain with trocar					
SB500-51	3x390x800mm	10FG			Box 10
SB500-52	5x390x800mm	15FG			Box 10
SB500-53	7x390x800mm	19FG			Box 10
SB500-53ST	small trocar 7x390x800mm	19FG			Box 10
SB500-54	8x390x800mm	flat			Box 10
SB300-61	round perforated 4x390x900mm	12FG			Box 10
Surgiblake silicone wound drain without trocar					
SB500-41	3x390x800mm	10FG			Box 10
SB500-42	5x390x800mm	15FG			Box 10
SB500-43	11x390x800mm	19FG			Box 10
SB500-46	8x350x800mm	24FG			Box 10
SB500-44	8x390x800mm	flat			Box 10
SB500-45	11x390x800mm	flat			Box 10
SB500-37	hubless flat JP style	7mm			Box 10
SB500-38	hubless flat JP style	11mm			Box 10
Surgiblake bulb drains					
SB100-02	100ml bulb, trocar	10FG			Box 10
SB100-12	100ml bulb only	-	-		Box 10
SB500-80	100ml bulb	flat perforated 10mm			Box 10
SB500-81	400ml bulb	flat perforated 10mm			Box 10
SB100-40	400ml bulb only	-	-		Box 10
Bendable trocars					
SB200-10	-	10FG	-		Box 10
SB200-15	-	15FG	-		Box 10

## CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au

Contact your National Surgical Corporation

Account Manager for more information



BE SUPPORTED BY AN  
AUSTRALIAN COMPANY  
PROVIDING VALUE BEYOND THE PRODUCT

**NATIONAL SURGICAL CORPORATION**

A proudly founded and owned Australian company - providing value to our customers beyond the product