

Medichill Instant Ice Packs

Recommended for use in all Emergency Departments and wards without refrigeration facilities as NO REFRIGERATION REQUIRED.

Three simple steps to activate, 1. Punch or squeeze pack 2. Shake 3. Apply. The ice pack will drop in temperature to around 0 degrees, which is dependent on the ambient temperatures at the time. Two sizes available, Regular and Large. No shelf life restriction, they will last for years and will always therefore be ready for use when needed most. Easy to use and suitable for adults or children, Pack has rounded corners for safety. AMMONIUM NITRATE FREE. TGA Registered 156748.

Covers can be purchased separately for additional comfort Product Code IPC32/42-100.



A carton contains 50 Ice Packs

**Contact your NSC
Account Manager
for product samples**




NATIONAL SURGICAL CORPORATION

*A proudly Australian founded and owned company:
providing value to the customer beyond the product*



Easy to use and apply



Contact your NSC Account Manager for more information or product samples



National Surgical Corporation Customer Service

1800 138 138 info@nationalsurgical.com.au

Visit us at www.nationalsurgical.com.au

