



**NATIONAL SURGICAL CORPORATION**

# CSSD Cleaning Brushes & Sterilisation Accessories

National Surgical Corporation (NSC) have a expert understanding of the requirements of CSSD departments and provide a range of essential products to improve efficiency in the cleaning and sterilisation process.

Our range of CSSD cleaning brushes and accessories have been selected from the highest quality manufacturers and our entire range is available for your CSSD department.

Be supported by an Australian company - talk with NSC to see how we can support your CSSD department

*A proudly founded and owned Australian company - providing value to our customers beyond the product*

# Contents

<b>Product</b>	<b>Page #</b>	<b>Product</b>	<b>Page #</b>
Book - sterilisation guide/textbook	24	Multi use cleaning brushes	4
Bottle brushes	12	Nylon twisted wire cleaning brushes	5, 6
Cannula cleaning brushes	6, 8- 11	Osteotome guards/ protectors	16, 17
Channel cleaning brushes	6, 8- 11	Peach clamp/clip	13
Coaching cards	24	Reamer cleaning brushes	7
Endoscopy cleaning brushes	8- 11	Scope cleaning brushes .	8- 11
Eye shields	29	Sharpness test material	26
Face shields	28, 29	Silicone mats	22
Flat cleaning brushes	13	Stainless steel bristle cleaning brushes	4
Flexible cleaning brushes	5 - 11	Steriliser cleaning set	12
General use cleaning brushes	4	Sterilisation trays and baskets	20
Hand and nail brushes	26	Stringers	18
Inspection mat	25	Suction tube cleaning brushes	6
Instrument tip protectors	16, 17	Surgical instrument/ stain remover pad	7
Jar brushes	12	Training and education	24
Keyboards	27	Tray corner guards	23
Kits/ set	6	Trays contamination	14
Large instrument brush	13	Valve cleaning brushes	11
Locking tags/Communication tags	19	Wall mounts	5

Face & eye shields




PAGE  
**28**

Scope cleaning brushes



PAGE  
**8**

Stringers




PAGE  
**18**

Communication/locking tags



PAGE  
**19**

Nylon twisted wire brushes



PAGE  
**5**

Steriliser cleaning kit



PAGE  
**12**

# Instrument Cleaning Brushes

## Cleaning brush selection guide

Having the correct cleaning brush is essential for providing a properly cleaned and sterilised instrument, ready for use again in the operating theatre. Most instruments will come into the CSSD with a level of blood, debris and other matter arising from surgery, it is vital that this is thoroughly cleaned and inspected prior to sterilisation.

*“If it’s not clean, it can’t be effectively sterilised”*

The type of brush that will be used will depend on the instrument being cleaned. Different brushes are available to provide the best option for cleaning the instrument.

## Toothbrush style cleaning brushes

The #1 brush used by all sterilisation departments, this brush comes in single or double ended. Nylon bristles are most commonly used, however, stainless steel and brass bristles are also available. These brushes are used to clean the external surfaces and box locks of scissors, forceps, needle holders and retractors. See Page 4



## Scope/channel cleaning brushes

Channel cleaning brushes are used for cleaning lumened instruments such as cannulas, laparoscopic instruments, scopes and suction tubes. These are difficult to clean as they have a narrow or small channel and have the largest risk of not being cleaned effectively.

### The “2mm” rule

The diameter of the brush should be 2mm larger than the cannula eg. 5mm cannula should be used with a 7mm bristle head brush. This allows for correct contact with the surface of the channel. If it’s too small, it won’t make contact.

Scope brushes are available in long lengths, with single or double ended bristles and are used as pull through brushes. The bristle head should always exit the cannula end completely, each time the brush is pushed into the channel/cannula.

Brushes to use: Nylon twisted wire brushes, blue/smooth tip brushes, channel brushes



## Cleaning brush guide

Instrument	Brush recommended	Page #
Baron & frazier suction tubes	Nylon bristles, wire shaft, flexible brushes	6, 9, 10, 11
Closed or semi-closed instruments	Fan tipped brushes	10
Instrument with short channels	Tapered bristles brushes	10
External instruments	Toothbrush style brushes	4
Coated instruments	Nylon bristles, do not use stainless steel bristle brushes	4, 5, 7, 11, 13
Reamers	Valve and port double ended brushes	7
Rasps	Circular head brushes	7
Laparoscopic instruments	Nylon bristles, twisted wire shafts	5, 6
Scopes	Nylon bristles, pull through brushes, long shaft lengths, flexible	5, 6, 7, 8, 9, 10, 11
Cannulated / channel instruments	Nylon bristles, pull through brushes, long shaft lengths, flexible	5, 6, 7, 8, 9, 10, 11

## Rust vs. Stains

Stains can be removed, whereas rust will leave permanent damage. To determine if a brown/orange discolouration is a stain or rust, use the eraser test. Rub a pencil eraser aggressively over the discolouration.

If the discolouration is removed with the eraser and the metal underneath is smooth and clean, this is a stain. If a pit mark appears under the discolouration, this corrosion is rust.





# Instrument Cleaning Brushes

## All Purpose, Multi-Use Cleaning Brushes

One end contains 3 rows of bristles for general use and the other end has a single row of bristles for delicate decontamination. These brushes are ideal for virtually cleaning any instrument; scissors, needle holders and forceps.

## Nylon Bristle Cleaning Brushes

For cleaning general instruments such as: scissors, hemostats, needle holders, dressing and tissue forceps. Plastic handle

Order #	Description	UOM
SH21650	Single ended nylon bristles 18cm	3/pk
SH54848	Double ended nylon bristles 18cm	3/pk
SH105606	Double ended nylon bristles 18cm firm, curved handle	3/pk
SH105607	Double ended nylon bristles 18cm firm, curved handle, half sized bristles	3/pk

CUSTOMER FAVOURITE



SH21650



SH54848



SH105606



SH105607

## Stainless Steel & Brass Bristle Cleaning Brushes

For cleaning general instruments such as: rasps, files, reamers, burrs, needle holders, jaws, vascular serrations, and stubborn stains. Plastic handle

Order #	Description	UOM
SH21651	Single ended stainless steel bristles 18cm	3/pk
SH54849	Double ended stainless steel bristles 18cm	3/pk
SH105604	Single ended brass bristles 18cm	3/pk



SH21651



SH54849



SH105604

## Double-Headed Cleaning Brush

Rigid, white, medical-grade nylon bristles and super-fine stainless steel bristles. Plastic handle

Order #	Description	UOM
SH105605	Double headed nylon & stainless steel bristles 18cm	3/pk



SH105605

# Instrument Cleaning Brushes



## Nylon Twisted Wire Cleaning Brushes

Can be used on: laparoscopic instruments, suction devices, and other general use instruments. Reusable

- White medical grade nylon bristles
- Made with anti-microbial nylon bristles and surgical stainless steel twisted and spun wire
- Designed in a way that could assist within the lumen during the cleaning process

CUSTOMER FAVOURITE

COLOUR INDICATES BRUSH DIAMETER



Order #	Description	UOM
<b>Nylon Twisted Wire Brush 30cm</b>		
SH45-542SS	2mm diameter, Red bead	3/pk
SH45-543SS	3mm diameter, White bead	3/pk
SH45-544-2SS	4mm diameter, Aqua bead	3/pk
SH45-545-2SS	5mm diameter, Green bead	3/pk
SH45-546-2SS	6mm diameter, Black bead	3/pk

Order #	Description	UOM
<b>Nylon Twisted Wire Brush 40.5cm</b>		
SH45-539SS	5mm diameter, Green bead	3/pk
SH45-540SS	7mm diameter, Yellow bead	3/pk
SH45-541SS	10mm diameter, Pink bead	3/pk
SH45-541-4SS	12mm diameter	3/pk
<b>Nylon Twisted Wire Brush 61cm</b>		
SH45-550SS	2.5mm diameter, Orange bead	3/pk
SH45-551SS	5mm diameter, Green bead	3/pk
SH45-552SS	7mm diameter, Yellow bead	3/pk
SH45-553SS	10mm diameter, Pink bead	3/pk



## Safely store your brushes with the wall mount rack!

- Reduce risk of damage to brushes
- Easily identify the diameter and length of the brush
- Can be custom labelled to suit the brushes stored
- Comes in a set of 2 racks
- Allows for easier reordering as codes are displayed

Order #	Description	UOM
RACK P	Brush wall rack/mount for storing nylon twisted wire brushes	2/set

## Continuous Bristles Multi Use Cleaning Brushes

Multi use cleaning brush - 2.5cm bristle interspaced every 5cm on 64cm long Stainless Steel wire





Order #	Description	UOM
<b>Delicate Scope Cleaning Brush 64cm</b>		
BR-IM-3	2.5mm diameter	3/pk
BR-IM-4	4mm diameter	3/pk
BR-IM-5	5mm diameter	3/pk

Order #	Description	UOM
<b>Delicate Scope Cleaning Brush</b>		
BR-IM-6	6mm diameter	3/pk
BR-IM-7	7mm diameter	3/pk
BR-IM-10	10mm diameter	3/pk

# Instrument Cleaning Brushes

## Smooth tip instrument cleaning brushes



With its long brush head and rigid design, these smooth ball tip nylon twisted brushes are ideal for cleaning many types of cannulated instruments.

- Long bristles provide more cleaning contact within the instrument
- Made with antimicrobial nylon bristles
- Some feature a blue ball tip on the end for added instrument protection
- Rigid design helps ensure ideal brush-cleaning capabilities
- Stainless-steel twisted wire and medical-grade nylon bristles

Order #	Description	Order #	Description
<b>Smooth tip nylon twisted wire brushes</b>		SH34317 ●	2.5mm diameter, 61cm
SH54844 ●	2mm diameter, 30.5cm	SH34318 ●	5mm diameter, 61cm
SH54845	3mm diameter, 30.5cm	SH34319	7mm diameter, 61cm
SH34306 ●	4mm diameter, 30.5cm	SH34320 ●	10mm diameter, 61cm
SH20659	5mm diameter, 30.5cm	SH34321	12mm diameter, 61cm
SH34308	6mm diameter, 30.5cm	<b>Smooth tip Frazier/Baron suction brushes</b>	
SH34309 ●	10mm diameter, 30.5cm	SH30618	3FG/ 1mm, 25cm
SH34310 ●	5mm diameter, 40.6cm	SH54759 ●	5FG/ 1.67mm, 25cm
SH34311 ●	6mm diameter, 40.6cm	SH30620	6FG/ 2mm, 25cm
SH34312 ●	7mm diameter, 40.6cm	SH30621 ●	7FG/ 2.3mm, 25cm
SH34313 ●	10mm diameter, 40.6cm	SH30622 ●	8FG, 2.7mm 25cm
SH34314 ●	12mm diameter, 40.6cm	SH30623	9FG/ 3mm, 25cm
SH34315 ●	15mm diameter, 40.6cm	SH30624	10FG/ 3.3mm, 25cm
SH34316 ●	20mm diameter, 40.6cm	SH30625	11FG/ 3.7mm, 25cm
		SH30626 ●	12FG/ 4mm, 25cm

### Essential brush set

Code: SH103631

Includes 30 brushes

1 pack (3) of each of the following brushes:

Order #	Description
<b>Smooth tip cleaning brushes</b>	
SH54845	3mm diameter, 30.5cm
SH34306	4mm diameter, 30.5cm
SH34308	6mm diameter, 30.5cm
SH34310	5mm diameter, 40.6cm
SH34312	7mm diameter, 40.6cm
SH34313	10mm diameter, 40.6cm
SH34318	5mm diameter, 61cm
<b>Frazier suction cleaning brush</b>	
SH54759	5FG, 25cm
<b>Double ended toothbrush style cleaning brush</b>	
SH54848	Nylon bristle
SH54849	Stainless steel bristle

### Starter set

Code: SH103630

Includes 39 brushes

1 pack (3) of each of the following brushes:

Order #	Description
<b>Smooth tip cleaning brushes</b>	
SH54844	2mm diameter, 30.5cm
SH54845	3mm diameter, 30.5cm
SH34306	4mm diameter, 30.5cm
SH20659	5mm diameter, 30.5cm
SH34308	6mm diameter, 30.5cm
SH34310	5mm diameter, 40.6cm
SH34312	7mm diameter, 40.6cm
SH34313	10mm diameter, 40.6cm
SH34314	12mm diameter, 61cm
SH34317	2.5mm diameter, 61cm
SH34318	5mm diameter, 61cm
SH34319	7mm diameter, 61cm
SH34320	10mm diameter, 61cm

# Instrument Cleaning Brushes

## For Hard to Clean Instruments



- Extra stiff medical-grade anti-microbial nylon bristles
- Aggressive cleaning
- Thicker, more rigid shafts to prevent bending
- Surgical grade stainless steel shafts

Order #	Description	UOM
<b>Blue nylon bristles brushes 46cm</b>		
BR-45-447	5mm diameter	3/pk
BR-45-448	10mm diameter	3/pk
BR-45-449	15mm diameter	3/pk
<b>White nylon bristles brushes 46cm</b>		
BR-45-667	5mm diameter	3/pk
BR-45-668	10mm diameter	3/pk
BR-45-669	15mm diameter	3/pk

## Stiff-bristled brush cleans cement and bone particles



Excellent for cleaning orthopaedic rasps and other difficult orthopaedic instruments, the head of this cleaning brush is designed so the rasp can pass through it, cleaning both sides at the same time.

- Antibacterial polyamide bristles which prevent or slow bacterial growth between two cleanings or brush sterilisations
- Ergonomic, steamable thermoplastic handle
- Reusable
- Long-lasting design withstands numerous sterilisation cycles
- Latex free
- Size: 23cm long overall; brush head is 2x7x9cm

Order #	Description	UOM
SH103697	Rasps orthopaedic cleaning brush 23cm	3/pk

## Surgical instrument stain/rust remover pad

This pad is the professional choice for removing all stains on surgical instruments.

The pads readily remove surface contamination and are used by the leading surgical instrument manufacturing and repair companies throughout the world

Order #	Description	UOM
SH22175	Stain remover pad, 15cm x 11cm	3/pk



Additional use:  
cleaning needles,  
ophthalmology

Caution: Do not  
use on any coated  
instruments

## Reamer antimicrobial cleaning brush

This double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers.

The circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Durable and reusable.

Specifications:

- Circular brush head diameter = 45mm, length = 50mm
- Straight brush end diameter = 8mm tapered, length = 30mm
- Overall length = 16.5cm


Order #	Description	UOM
SH45-755SS	Reamer cleaning brush double ended 16.5cm	3/pk



# Scope/Channel Cleaning Brushes

## Disposable, Flexible Channel Brushes

A full range of disposable brushes ideal for cleaning flexible endoscopy equipment; including but not limited to bronchoscopes, colonoscopes, duodenoscopes, gastroscopes, rhinoscopes and sigmoidoscopes and other general use instruments.



Order #	Description	UOM
<b>75cm Brushes</b>		
SH18898	1.5mm, Blue, 4.5FG	15/pk
SH18897	2.4mm, White, 7FG	15/pk
SH18896	3.2mm, Orange, 10FG	15/pk
<b>100cm Brushes</b>		
SH21479	1.5mm, Blue, 4.5FG	15/pk
SH21478	2.4mm, White, 7FG	15/pk
SH21477	3.2mm, Orange, 10FG	15/pk
<b>150cm Brushes</b>		
SH21474	1.5mm, Blue, 4.5FG	15/pk
SH18902	2.4mm, White, 7FG	15/pk
SH21729	3.2mm, Orange, 9FG	15/pk
SH21731	4.2mm, Red, 13FG	15/pk
<b>200cm Brushes</b>		
SH18901	1.5mm, Blue, 4.5FG	15/pk
SH18900	2.4mm, White, 7FG	15/pk
SH21730	3.2mm, Orange, 9FG	15/pk
SH18899	4.2mm, Red, 13FG	15/pk

Call your  
NSC Account Manager for more  
information and pricing

**1800 138 138**

## External endoscope cleaning brushes

### Important brushes Ideal for an Important step in endoscope reprocessing



SH18441 - Reusable



SH18442 - Single use

These uniquely designed external cleaning brushes are ideal for cleaning the exterior of endoscopes and optical instruments, without the risk of damage.

With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments.

Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.

Brush length	Handle length	Total length
Single use - 70mm	8cm	15cm
Reusable - 90mm	12cm	21cm

Order #	Description	UOM
SH18441	Reusable - 12cm thermoplastic handle	3/pk
SH18442	Single use - 8cm handle with vinyl caps	50/pk



# Scope/Channel Cleaning Brushes



## Colour-coded sheath brushes for quick, easy identification

- Designed to preserve instrument integrity and last longer than twisted wire brushes
- Colour indicates brush diameter
- Rigid or flexible plastic sheath
- No sharp ends
- Won't scratch instruments
- Will not lose their intended shape

SHEATH COLOUR INDICATES BRUSH DIAMETER

2mm	3mm	4mm	5mm	6mm	8mm	10mm	12.5mm	15mm	20mm
-----	-----	-----	-----	-----	-----	------	--------	------	------


Description	Reusable Pack 3	Single Use Pack 20
ARC Plastic Sheath Brush - Flexible		
31.5cm x 2mm, White	SH62387	SH62427
33cm x 3mm, Yellow	SH62388	SH62428
33cm x 4mm, Violet	SH62389	SH62430
33cm x 5mm, Grey	SH62390	SH62431
46.5cm x 2mm, White	SH62395	SH62436
48cm x 3mm, Yellow	SH62396	SH62437
48cm x 4mm, Violet	SH62397	SH62439
48cm x 5mm, Grey	SH62398	SH62440
61cm x 2mm, White	SH62405	SH62447
61cm x 3mm, Yellow	SH62406	SH62448
61cm x 4mm, Violet	SH62407	SH62449
61cm x 5mm, Grey	SH62408	SH62450
76cm x 2mm, White	SH62415	SH62457
76cm x 3mm, Yellow	SH62416	SH62458
76cm x 4mm, Violet	SH62417	SH62459
76cm x 5mm, Grey	SH62418	SH62463

Description	Reusable Pack 3	Single Use Pack 20
ARC Plastic Sheath Brush - Rigid		
35.5cm x 6mm, Green	SH62391	SH62432
35.5cm x 8mm, Blue	SH62392	SH62433
35.5cm x 10mm, Red	SH62393	SH62434
35.5cm x 12.5mm, Yellow	SH62394	SH62435
51cm x 6mm, Green	SH62399	SH62441
51cm x 8mm, Blue	SH62400	SH62442
51cm x 10mm, Red	SH62401	SH62443
51cm x 12.5mm, Yellow	SH62402	SH62444
51cm x 15mm, Green	SH62403	SH62445
51cm x 20mm, Blue	SH62404	SH62446
61cm x 6mm, Green	SH62409	SH62451
61cm x 8mm, Blue	SH62410	SH62452
61cm x 10mm, Red	SH62411	SH62453
61cm x 12.5mm, Yellow	SH62412	SH62454
61cm x 15mm, Green	SH62413	SH62455
61cm x 20mm, Blue	SH62414	SH62456
76cm x 5mm, Rigid White	SH62419	SH62464
76cm x 6mm, Green	SH62420	SH62465
76cm x 8mm, Blue	SH62421	SH62466
76cm x 10mm, Red	SH62422	SH62467
80cm x 12.5mm, Yellow	SH62424	SH62468
80cm x 15mm, Green	SH62425	SH62469
80cm x 20mm, Blue	SH62426	SH62470

Smaller diameter brushes have more of a flexible sheath for help navigating fine channels.

A more rigid sheath is used with larger diameter brushes to provide more pushing power

# Scope/Channel Cleaning Brushes




## DayBrush™ cylindrical channel cleaning

- Bristle fan tip protects and cleans channels
- Rigid Delrin is easy to insert and manipulate, helping to reduce time
- No sharp ends or hard metal tips that can scratch or cause damage
- Ideal for cleaning endoscopes (up to 10 scopes before discarding)
- Discard after cleaning 10 scopes or sooner if any brush degradation is noticed
- Sheath = 1mm, Brush dia = 2mm, Brush length = 15mm

Order #	Description	UOM
SH22173	Single ended channel brush 230cm x 2mm dia	50/pk

## Endoscopy channel cleaning brushes - tapered and cylindrical



- Long brush head (15-24mm) has more bristles for more thorough cleaning within the channel
- Bristle fan tip protects and cleans the endoscope channel
- No sharp ends of hard metal tips that can scratch or cause damage
- Sheaths: dia 3mm = 1.55mm, dia 5mm = 1.75mm, tapered = 1.85mm

Order #	Description	UOM
<b>Single ended channel brush</b>		
SH22172	Delrin 100cm x 5mm dia, rigid, metallic grey	50/pk
SH22167	Polypropylene, 100cm x 5mm dia, flexible, pink	50/pk
SH22171	180cm x 3mm dia, purple	50/pk
SH22166	180cm x 5mm dia, pink	50/pk
SH22170	230cm x 3mm dia, purple	50/pk
SH22162	230cm x 5mm dia, pink	50/pk
SH22164	230cm x 5mm dia incl valve brush double ended, pink	50/pk
SH22168	Tapered bristles 5-6mm x 230cm, pink	50/pk
<b>Double ended channel brush</b>		
SH22163	230cm x 5mm dia, pink	50/pk
SH22165	230cm x 5mm dia incl valve brush double ended, red	50/pk
SH22169	Tapered bristles 5-6mm x 230cm, pink	50/pk

# Scope/Channel Cleaning Brushes

## Precision external brush

Precise and quick to remove difficult debris, the reusable precision brush is a compact external cleaning brush that eliminates the need for the stainless-steel toothbrush style brush.

The result is reduced risk of instrument damage and personal injury.

- Great for easily cleaning deep into instrument jaws
- Superior cleaning strength with polyamide TYNEX® bristles
- Steamable, thermoplastic handle
- Extra-firm bristle density
- Steam/chemical sterilisable
- Handle length = 120mm, brush length = 40mm  
brush width = 12mm, brush thickness = 8mm




Order #	Description	UOM
SH18443	Reusable precision brush 16cm	3/pk

## Double-ended trumpet valve brushes

- Higher density of bristles translates into more cleaning contact within the channel
- Bristle fan tip on the head protects and cleans the channel
- No sharp ends or hard metal tips that can scratch or cause damage
- Can be used to clean ports and valves of endoscopes

Order #	Description	UOM
SH18621	White 5/7mm, 12.5cm	3/pk
SH27029	White 5/7mm, 12.5cm	50/pk
SH27032	Yellow 5/11mm, 15cm	3/pk
SH27033	Yellow 5/11mm, 15cm	50/pk





## Endoscope cleaning care kit

Convenient endoscope cleaning kit, includes brushes for cleaning and a sponge for wiping. Ideal for use with most popular flexible endoscopes

Kit includes:

- Valve brush double ended 5/11mm, 15cm long (SH27032)
- Channel cleaning brush 230cm long x either 2mm or 3mm diameter
- Safe-Soak endo sponge

Order #	Description	UOM
SH105602	Set with 3mm diameter brush (purple)	40/sets
SH105603	Set with 2mm diameter brush (yellow)	40/sets



# General Use Cleaning Brushes

## Autoclave cleaning brush kit

Ideal for cleaning autoclaves/sterilisers, and other general use instruments. Medical grade nylon bristles

- Improves steriliser operating efficiency
- Improves drying efficiency

**"Everything you need  
to maintain your  
steriliser"**



Order #	Description	UOM
SH34088	Autoclave/Steriliser Cleaning Brush Set	Set 3
<b>Replacement parts</b>		
SH34087	Brush head for use with autoclave extension pole	1/ea
SH34089	Anodised aluminium extension pole 60 - 122cm	1/ea

### SH34088

#### Cleaning Brush Set Includes:

- 1 x Handheld Utility Brush 9x23cm
- 1 x Brush Handle/ extension pole 60-122cm
- 1 x Rounded Corner Brush Head 14x20cm

## Jar/bottle brush



- Long, flexible bristles to maximise large surface cleaning and ensure a wider array of cleaning applications
- Lightweight ergonomic handle
- Autoclavable
- Designed specifically for hospital applications
- High-quality polyamide clear bristles make it easier to spot debris

Order #	Description	UOM
SH20056	Reusable External Cleaning Brush 35.5cm, brush head - 7cm dia x 12cm	3/pk

## Large head cleaning brush with handle

Utility brush can be used on: sinks, autoclaves, ultrasonic cleaners, bowls, basins, I.V. poles, commodes, wheelchairs, and other general use products.

9cm x 23cm including sturdy plastic handle

Order #	Description	UOM
SH34086	Large head cleaning brush 9x23cm	1/ea





## Large instrument cleaning brush



Can be used on: retractors, bone cutters, specula, rib spreaders, bowls, basins and other general use instruments.

- Plastic handle
- Nylon bristles
- Overall length 26cm (2cm hanging hook on end)
- Handle length 12cm

Order #	Description	UOM
694B	Plastic handle nylon brush 26cm	3/pk

## Theatre/peach clip, clamp

Heavy-duty clamp is designed to secure drapes to drape stands, IV poles, etc. Can be used for a variety of clamping needs.

- Washable (not sterilisable), does not rust
- Size: 4.75cm wide x 7cm long
- Comes in an assorted colours pack
- Non-sterile



Order #	Description	UOM
<b>Peach Clip/ Clamp</b>		
PP-JPC-3NS	Standard pack, assorted colours	8/pk
PP-JPC-3NSB	Big pack, assorted colours	16/pk
PP-JPC-1086-B50	Big pack, blue only	50/pk

## Flat bristle brushes

- Long ergonomic handle for ease of use and strong grip
- Polyamide transparent fibre bristles
- Plastic handle will not release toxic components
- Available with soft or rigid bristles
- Extra long length of brush head 75mm



Order #	Description	UOM
<b>Reusable Flat Brush 22cm</b>		
SH23528	Blue - extra rigid bristles	3/pk
SH23489	Green - rigid bristles	3/pk
SH34084	Yellow - soft flexible bristles	3/pk
SH34085	Red - extra soft flexible bristles	3/pk

# Scope Protection & Accessories

## Probe & scope transport trays & cart

Safely and securely transport clean and dirty endoscopes and probes to and from OT and CSSD.

Comes with dividers to separate the probe from the plug and cable, while also protecting the delicate probe during transport.

- Included "Clean" & "Dirty" labels
- Clean with hospital-grade disinfectants and stackable
- Protect patients and personnel from cross-contamination
- Compatible with most endocavity probes, TEE probes, endoscopes

Nautilus transport tray storage cart

- Stores up to 5 Nautilus multi transport trays (code: SH43609)
- Lipped shelves prevent trays from sliding out
- Flat top is excellent for use as additional workspace
- Designed to move short distances or between rooms
- Made from lightweight and rigid, EPVC Sintra®, easy to clean

Order #	Description	UOM
SH43609	Nautilus-Multi transport tray with cover - 70x49x9cm	4/pk
SH41235	Probe transport tray with cover - 36x43x7.5cm	4/pk
SH87340	Nautilus transport tray storage cart - 89x53x94cm	1/ea



## Instrument tray system

- Designed to hold a complete instrument tray
- Helps reduce cross-contamination risks associated with touching individual instruments
- Snap-down lid and molded-in gasket provide a near watertight seal
- Max temperature: 94°C

Order #	Description	UOM
SH57587-RD	Instrument tray system - 36x63.5x15cm	1/ea

## Scope lens tip protectors

Protect endoscope tips and lenses during transport and storage

- Allows moisture to escape eliminating potential bacteria growth
- Semi-rigid expandable mesh protects the tip/lens
- Securely fits many different sizes of scopes and probes
- Protectors can stay on the scope tip until used and then put back on after use for transport back



Order #	Description	UOM
<b>Tip Protectors for Scope Lens</b>		
SH43884	Small Red, 2 - 10mm	100/pk
SH43885	Large Blue, 7 - 15mm	100/pk

## High resistance trays (yellow)


- Protect your surgical instruments and equipment during transportation to/from Theatre and CSSD
- Machine washable and sterilisable
- Increased chemical resistance to most common hospital detergents
- Can be used with sealed tag to mark contaminated instruments during transportation

Order #	Description	UOM
T104060-ND-HR	Tray high resistance - 60x40x10cm	1/ea
L4060-HR	Lid for sterile high resistance tray	1/ea

## Anti-microbial trays (green)

- Green infection control trays can be used in areas where patients are at enhanced risk for contracting infections
- These trays resist the growth of microbial and mould throughout the entire lifespan of the module
- Proven to reduce the levels of microbial on the module surface by up to 99% within 24 hours
- Max temperature: 95°C

Order #	Description	UOM
T104060-AM	Tray antimicrobial - 60x40x10cm	1/ea
T204060-AM	Tray antimicrobial - 60x40x20cm	1/ea



High resistance with lid



T104060-ND-HR + LID



T204060-AM



T104060-AM

## EndoSock™ flexible scope sleeve



- Protect scopes from bacteria between uses
- Paper-based sleeve protects against moisture
- Creates a barrier against physical contact
- Open on both ends for easy application and removal

Order #	Description	UOM
<b>EndoSock™ Flexible Scope Sleeve</b>		
SH41570	Small 91cm	200/pk
SH41571	Large 198cm	100/pk

## Endoscope tip guards, eye piece

Order #	Description	UOM
SH20271	Endo tip guards, non-vented	6/pk



# Instrument Tip Protection

Instrument tip protectors are designed to protect medical instruments from damage during sterilisation, storage, and handling



## Round cap guards - vented, clear

Order #	Description	UOM
SH20195	1.6 x 19mm	100/pk
SH20196	2 x 19mm	100/pk
SH20197	2.8 x 19mm	100/pk
SH20198	3.2 x 25.4mm	100/pk
SH20199	5 x 25mm	100/pk
SH20200	10 x 19mm	100/pk
SH20201	2 x 9 x 25mm	100/pk
SH20202	2 x 16 x 25mm	100/pk
SH20203	3 x 25 x 25mm	100/pk
SH20204	Assorted	100/pk



## Weitlaner retractor protectors

Order #	Description	UOM
SH54749	Small pink	50/pk
SH54753	Small clear	50/pk
SH54754	Medium green	50/pk
SH54755	Medium clear	50/pk

## Vented laparoscopic tip protectors

Order #	Description	UOM
SH20298	20 x 5mm	25/pk
SH20299	43 x 5mm	20/pk
SH20300	20 x 10mm	25/pk
SH20301	57 x 10mm • Trocar probe	20/pk
SH20297	19mm • Verres needle protector	50/pk



## Retractor tip guards

Vented, clear for protection of retractor tips



Order #	Description	UOM
SH54756	38mm	55/pk
SH54757	51mm	45/pk
SH54758	76mm	45/pk

## Instrument tip protectors

Order #	Description	UOM	Order #	Description	UOM
<b>Vented Instrument Guards</b>			<b>Non-Vented Instrument Guards</b>		
SH54692	1.6 x 19mm • White	100/pk	SH54682	1.6 x 19mm • White	100/pk
SH54693	2 x 19mm • Blue	100/pk	SH54683	2 x 19mm • Blue	100/pk
SH54694	2.8 x 19mm • Green	100/pk	SH54684	2.8 x 19mm • Green	100/pk
SH54695	3.2 x 25.4mm • Maroon	100/pk	SH54685	3.2 x 25.4mm • Maroon	100/pk
SH54696	4.8 x 25.4mm • Yellow	100/pk	SH54686	4.8 x 25.4mm • Yellow	100/pk
SH54697	10.3 x 19mm • Purple	100/pk	SH54687	10.3 x 19mm • Purple	100/pk
SH54698	1.6 x 9.5 x 25.4mm • Orange	100/pk	SH54688	1.6 x 9.5 x 25.4mm • Orange	100/pk
SH54699	2 x 15.9 x 25.4mm • Brown	100/pk	SH54689	2 x 15.9 x 25.4mm • Brown	100/pk
SH54702	3.2 x 25.4 x 25.4mm • Black	100/pk	SH54690	3.2 x 25.4 x 25.4mm • Black	100/pk



All tip protectors on these pages are compatible with steam sterilisation and conventional hospital ethylene oxide. Single use, radio-opaque, latex free



# Instrument Tip Protection

## Webbed osteotome guards



Order #	Description	UOM
SH20302	10 x 76mm, Orange	24/pk
SH20303	15 x 76mm, Green	24/pk
SH20304	20 x 76mm, Blue	24/pk
SH20305	40 x 76mm, Red	24/pk
SH20306	60 x 76mm, Black	16/pk
SH20307	70 x 76mm, White	16/pk

## Double perforated instrument tip protectors



Order #	Description	UOM
<b>Double perforated instrument tip protectors</b>		
SH54733	1.6 x 19mm • White	50/pk
SH54734	2 x 19mm • Blue	50/pk
SH54735	2.8 x 19mm • Green	50/pk
SH54736	3.2 x 25.4mm • Red	50/pk
SH54737	4.8 x 25.4mm • Yellow	50/pk
SH54738	1.6 x 9.5mm • Orange	50/pk
SH54739	Assorted	50/pk

## Standard osteotome tip guards



Description	Vented Order #	Non-Vented Order #	UOM
6.4 x 19mm	SH54740	SH20236	50/pk
9.5 x 25.4mm	SH54741	SH20237	50/pk
12.7 x 25.4mm	SH54742	SH20238	50/pk
19 x 25.4mm	SH54743	SH20239	50/pk
25 x 25.4mm	SH54744	SH20240	50/pk
32 x 32mm	SH54745	SH20241	50/pk
38 x 38mm	SH54746	SH20242	50/pk
44.5 x 32mm	SH54747	SH20243	50/pk
Assorted	SH54748		50/pk


## Clear tip guards



Description	Vented Order #	Non-Vented Order #	UOM
1.57 x 19mm	SH54718	SH54705	50/pk
2 x 19mm	SH54720	SH54706	50/pk
2.8 x 19mm	SH54721	SH54707	50/pk
3.2 x 25.4mm	SH54724	SH54708	50/pk
4.75 x 25.4mm	SH54725	SH54709	50/pk
10.3 x 19mm	SH54726	SH54712	50/pk
1.6 x 9.5 x 25.4mm	SH54727	SH54713	50/pk
2 x 16 x 25.4mm	SH54729	SH54714	50/pk
3.2 x 25.4 x 25.4mm	SH54731	SH54716	50/pk
Assorted	SH54732	SH54717	50/pk

# Instrument Stringers

Designed to keep hinged instruments in an open position during sterilisation. Stringers help keep instruments organised and maximise usable tray space. Stainless steel. Various options



Order #	Description	UOM
<b>A - Weinstein U-shaped stringers</b>		
SH111401	7.6x15.2cm	1/ea
SH111402	7.6x20.3cm	1/ea
SH111403	7.6x25.4cm	1/ea
SH111404	7.6x30.5cm	1/ea
SH111405	7.6x35.6cm	1/ea
SH111406	7.6x40.6cm	1/ea
<b>B - Saunders</b>		
SH111408	Spring wire rack	1/ea
<b>C - Bunt forceps holder with catch</b>		
SH111409	12.1cm	6/pk
<b>D - "T" Bar instrument holders 15.2cm</b>		
SH111410	15.2cm	1/ea
SH111411	30.5cm	1/ea
<b>E - Pin-Lok instrument stringers</b>		
SH111412	6.4x10.2cm	1/ea
SH111413	6.4x15.2cm	1/ea
SH111414	6.4x20.3cm	1/ea
SH111415	6.4x25.4cm	1/ea
SH111416	6.4x30.5cm	1/ea
SH111417	6.4x35.6cm	1/ea
SH111418	6.4x40.6cm	1/ea
SH111419	6.4x45.7cm	1/ea
SH111420	6.4x50.8cm	1/ea
SH111421	8.9x45.7cm	1/ea
SH111422	8.9x50.8cm	1/ea

Order #	Description	UOM
<b>F - Open top stringers</b>		
SH111423	6.4x15.2cm	1/ea
SH111424	6.4x20.3cm	1/ea
SH111425	6.4x25.4cm	1/ea
SH111426	6.4x30.5cm	1/ea
SH111427	6.4x35.6cm	1/ea
SH111428	6.4x40.6cm	1/ea
SH111429	6.4x45.7cm	1/ea
<b>G - Hook-lok stringers</b>		
SH111430	6.4x15.2cm	1/ea
SH111431	6.4x20.3cm	1/ea
SH111432	6.4x25.4cm	1/ea
SH111433	6.4x30.5cm	1/ea
SH111434	6.4x35.6cm	1/ea
SH111435	6.4x40.6cm	1/ea
<b>H - Side opening</b>		
ACC-5010	10x6.4cm	1/ea
ACC-5011	15.2x6.4cm	1/ea
ACC-5012	20.3x6.4cm	1/ea
ACC-5013	25.4x6.4cm	1/ea
ACC-5014	30.5x6.4cm	1/ea

Order #	Description	UOM
<b>I - Centre hinge lok stringers</b>		
SH111437	6.4x15.2cm	1/ea
SH111439	6.4x25.4cm	1/ea
SH111440	6.4x30.5cm	1/ea
SH111441	6.4x35.6cm	1/ea
SH111442	6.4x40.6cm	1/ea
SH111443	6.4x45.7cm	1/ea
<b>J - Centre Male/Female lok stringers</b>		
SH111444	6.4x15.2cm	1/ea
SH111445	6.4x20.3cm	1/ea
SH111446	6.4x25.4cm	1/ea
SH111447	6.4x30.5cm	1/ea
SH111448	6.4x35.6cm	1/ea
SH111449	6.4x40.6cm	1/ea
SH111450	6.4x45.7cm	1/ea

## Facilitate faster, safer instrument handling

- Unique design promotes faster, safer instrument handling
- Allows for fast two-handed instrument exchange
- Use with instrument stringers of any width
- High-density polyethylene
- PETG base enhances stability
- Bottom bumpers help keep the stand in place

Order #	Description	UOM
SH79494	Instrument stringer stand	1/ea



# Communication, locking tags

## Improve communication between Theatre and CSSD



Easytwist indicator dot



Easytwist solid colour



Easytwist printed



square style printed



square style solid colour

### Solid colour • Pre-printed label • Sterilisation indicator dot

- Use on sterilisation containers, trays, individual instruments, carts or to secure goods
- Prevent tray, equipment, and instrument tampering
- Multiple colour options are ideal for colour-coding instruments and equipment throughout your facility
- Large base for ease of handling and application
- EasyTwist breaking or pull out design
- Fully open (45°) to facilitate insertion into instrument containers
- Will not melt or soften with any steam-sterilisation protocols
- DEHP-free polypropylene



Order #	Description	UOM	Order #	Description	UOM
<b>Locking tag solid colour square style</b>			<b>Locking tag indicator dot Easytwist</b>		
LT-100	White	100/pk	SH42407-BK	Black	100/pk
LT-110	Pink	100/pk	SH42407-BL	Blue	100/pk
LT-112	Gold	100/pk	SH42407-GN	Green	100/pk
LT-116	Purple	100/pk	SH42407-OR	Orange	100/pk
LT-200	Blue	100/pk	SH42407-PK	Pink	100/pk
LT-300	Green	100/pk	SH42407-RD	Red	100/pk
LT-400	Red	100/pk	SH42407-WH	White	100/pk
LT-500	Yellow	100/pk	<b>Locking tag printed Easytwist</b>		
LT-600	Grey	100/pk	SH42408	"Clean" White	100/pk
LT-700	Orange	100/pk	SH42409	"Sterile" Blue	100/pk
LT-800	Brown	100/pk	SH42410	"Repair" Red	100/pk
LT-900	Black	100/pk	SH42411	"Sharpen" Yellow	100/pk
<b>Locking tag solid colour Easytwist</b>			SH42412	"Loaner" Yellow	100/pk
SH42406-BK	Black	100/pk	SH42413	"Discard" Black	100/pk
SH42406-BL	Blue	100/pk	<b>Locking tag printed square style</b>		
SH42406-GN	Green	100/pk	LT-210	"Quick Turnover" Blue	100/pk
SH42406-OR	Orange	100/pk	LT-310	"Loaner" Green	100/pk
SH42406-PK	Pink	100/pk	LT-450	"Repair" Red	100/pk
SH42406-RD	Red	100/pk	LT-550	"Sharpen" Yellow	100/pk
SH42406-WH	White	100/pk	LT-720	"Broken" Orange	100/pk
			LT-930	"Discard" Black	100/pk

# Sterilisation Trays & Baskets

## Stainless steel wire sterilisation baskets

NSC sterilisation trays and baskets are high quality, manufactured in electrolytically polished 304 stainless steel. They can be used to securely store and transport different sized instruments for cleaning and sterilisation. DIN sizes available. NSC can custom make any size.


See all the sizes and styles available in the NSC Sterilisation Baskets & Trays brochure

### Perforated base, wire mesh sides

PWB-407  
240x250x50mm



PWB-451  
480x250x50mm



PWB-470  
300x270x60mm



PWB-469  
270x150x30mm



PWB-492  
270x150x40mm

PWB-449  
540x250x70mm



PWB-423  
405x250x60mm



PWB-453  
480x250x100mm



PWB-448  
540x250x60mm



### Woven wire mesh base & sides

WMB-607  
240x250x50mm



WMB-651  
480x250x50mm



WMB-653  
480x250x100mm



Lid removable for wire mesh and perforated baskets - various sizes



### Micro mesh baskets with lids 2x2mm mesh

MMB-729  
40x40x20mm



MMB-731  
80x40x20mm



### Super fine mesh baskets with lids 1x1mm mesh

SMB-810  
40x40x20mm



SMB-815  
80x80x40mm



### Endoscope baskets with removable lids

ESB-903-CLIP  
290x80x52mm  
with clip/pin



ESB-905-CLIP  
460x80x52mm  
with clip/pin



### Instrument trays

Flat wire mesh base, perforated sides with silicone inserts and silicone strip





# Sterilisation Trays & Baskets

## Organiser/basket



Designed to organise and hold plastic/paper peel pouches during the sterilisation process

- Durable electro-polished stainless steel
- Provides spacing to ensure sufficient sterilant contact to all surfaces
- Provides drainage if moisture present
- Compatible with all sterilisation procedures
- Insert SW1-2000 fits easily into basket SWO-2000

Order #	Description	UOM
SW1-2000	Insert only for basket 36x51x11.5cm	1/ea
SWO-2000	Organiser basket (no insert) 36x53x16.5cm	1/ea
SW2-2000	Basket/Insert Combo 36x53x16.5cm	1/ea

## Kirschner wire Stainless Steel sterilisation tube

Tube container 18/8 stainless steel for holding Kirschner wires up to 15cm long.

Compatible with all sterilisation procedures. Provides excellent drainage in tube. Stainless steel cap provides protection for wire tips and trocar ends.

Container size: dia. 12mm




*supplied without wires*

Order #	Description	UOM
<b>Kirschner wire tube Ø12mm</b>		
76.16.03	16cm for wires up to 15cm long	1/ea
76.16.04	19cm for wires up to 18cm long	1/ea
76.16.05	32cm for wires up to 31cm long	1/ea

## Sterilisation racks



CSDD-SIR-1



CSDD-SIR-2

Ideal to securely hold pouched instruments during the sterilisation process.

SIR-2 comes in a pair - one rack sits inside the other so that rack width can be adjusted wider by sliding sideways


Order #	Description	UOM
CSDD-SIR-1	Steriliser rack - 130x80x50mm	1/ea
CSDD-SIR-2	Steriliser rack - 890x150x200mm	2/pc

## Silicone inserts for sterilisation trays

Silicone inserts/bars can be manufactured in the tray when ordered or available separately and attached with screws to any tray.

Can be moved to any position on tray using screws and secure base or available with push through base

Custom options available



holes in nibs



solid nibs



square with holes



different sizes available



push through base options



removable base



round style



scope holders


# Sterilisation Tray Accessories

## Silicone mats

Can be used in any sterilisation tray or basket. Can be cut down to any size. Perforated to allow exposure and ventilation during sterilisation. Different square hole options

- Thermo-resistant from - 60°C to + 300°C
- 100% sterilisable and washable

Order #	Description	UOM
<b>Silicone mat - holes spaced 25mm, hole Ø 5mm</b>		
ACC-5015	515x240x15mm	1/ea
ACC-5015-SP	250x240x15mm	1/ea
ACC-5016	545x565x19mm	1/ea
<b>Silicone mat - 330x460mm</b>		
SH18444	Single pack	1/ea
SH18445	Bulk pack	5/pk
<b>Vicolab® silicone mats</b>		
SHV558908	500x320mm	1/ea
SHV569104	For micro instruments - 500x320mm	1/ea



## Silicone flat mat - use as a cover or base mat

Protects your instruments and secures instruments in baskets from moving around during wash cycles. Can also be placed under the instruments

Specially designed mats for trays used in autoclaves and ultrasonic cleaners. 15mm diameter large square holes to allow for water to flow out

- Sterilisable, washable
- Fast dry
- Thermo-resistant from - 60°C to + 300°C




Order #	Description	UOM
SHV588690	Vicolab® silicone flat mat - 540 x320mm	1/ea

## Non-slip silicone mat

Place under stored sterile trays and stop them from moving and potentially contaminating sterile instruments.

Made from Silicone (HCR) High Consistency Silicone Rubber, can be sterilised, place cleaned trays on the silicone mat and sterilise with the tray. Lifespan of 100,000 cycles



Order #	Description	UOM
SIMAT-250X420	420x250mm	1/ea
SIMAT-250X255	255x250mm	1/ea

## Sterilisation tray corner guards

- Prevents tearing caused by sharp corners or feet on surgical trays
- Fits virtually all types of tray corners, including rounded ones
- Attention-grabbing bright-orange colour makes them easy to identify and helps reduce their risk of being discarded
- Withstands repeated sterilisation up to 134°C
- Reduces the need to reprocess surgical trays due to ripped sterile wrapping
- Helps minimise time and costs associated with reprocessing sterile surgical trays
- Elevate height gives your trays space to help circulate air



SH54677



SH54676



SH54675



CP-01/SH16975


Order #	Description	UOM
SH54677	SoftPockets™ foam corner protectors	80/Sets of 4
SH54676	Corner protectors disposable blue	24/pk
SH54675	Corner protectors reusable orange/red	24/pk
CP-01	Corner protectors silicone pink	8/pk
SH16975		80/pk

## Sterilisable work surface mat

The flat mat is ideal for protection of surgical instruments whilst cleaning and inspecting them. Working on a stainless steel surface can be hard on the instruments. With the silicone mat you can protect sensitive pointed instruments and preserve sharp, fragile and delicate instruments against continually aggressive knocks, increasing the length of the lifetime of the instruments.

- Noise is diminished because there is no contact of instruments direct on the stainless steel worktable
- No reflective glare off the stainless steel work table due to coverage of silicone flat mat, working conditions are improved
- Large size: 2400x1200x3mm
- Cut to any size, for example to put around microscopes areas
- Thermo-resistant (- 60° to + 300°C)

Order #	Description	UOM
SHV578623	Silicone work surface mat - 2400x1200x3mm	1/ea





# CSSD Personnel Education

## Instrument coaching cards for training and continuing education



- Full-colour, detailed photographs on every card
- Includes proper instrument name and important points of inspection
- Learn sharpness test standards
- 163 double-sided flashcards
- Helps increase technician's surgical tray accuracy

Order #	Description	UOM
SH90861	Instrument coaching cards	163/pk

## "The World of Surgical Instruments" Book

### The Quintessential Universal Instrument Guide

- All-inclusive guide includes surgical-instrument names, inspection points, sharpness testing, and measuring
- Excellent study guide and resource for healthcare professionals, CSSD personnel, technicians and managers
- Increases patient safety and surgical satisfaction
- Over 1,100 high-resolution images
- Step-by-step approach for ease of comprehension

Order #	Description	UOM
SH61826	The World of Surgical Instruments An Illustrated Guide Textbook	1/ea






# Sterilisation Department Essentials

## Surgical instrument inspection mat

Heavy duty workplace mat features inspection tips for most common surgical instruments. Provides learning for new technicians and friendly reminders for seasoned professionals. Content provided by Rick Schultz, renowned author and instrument expert

- Learn proper names and inspection points of basic instruments
- Detailed diagrams assist with learning
- Durable construction
- Non-slip back
- Features 46cm ruler across the top, inch ruler along bottom

Order #	Description	UOM
SH105500	Instrument inspection mat - 66x43cm	1/ea



## 12-roll label tape dispenser

This label tape dispenser has a refined, innovative design which minimises hassles often experienced when removing several rolls at once to change out a single, empty roll. Unique four-arm “bays” hold up to three rolls on each arm.

- Two centre slots create four arms for easy exchanges
- Serrated blade for easy dispensing
- Made of clear polycarbonate with open areas for optimal cleaning
- Four suction cups secure placement on carts
- Mount on walls or slots using precut key holes
- Size: 33 x 8.7 x 5.5cm



Order #	Description	UOM
SH16603	Label tape dispenser for 12 rolls	1/ea

## Instrument lubricant spray bottle

Uniquely formulated with an All-Natural corrosion inhibitor. Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving.

Forms a protective coating which extends the instrument's life

- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

Order #	Description	UOM
SH104764	Instrument lubricant spray bottle 946mL	1/ea
SH105503		12/pk



# Sterilisation Department Essentials

## Scissor sharpness test material



- Quickly and easily test scissor sharpness
- Make several cuts in the material to test for sharpness
- Cut all the way to the distal tip when testing
- Convenient 46cm rolls
- Red testing material is designed for scissors 11cm or longer
- Yellow testing material is designed for scissors 11cm or shorter

Order #	Description	UOM
<b>Scissor Testing Material, 5.5m x 14cm Latex Free</b>		
SH39975	Yellow, Thin - for scissors 10.5cm or smaller	1/ea
SH39976	Red, Medium - for scissors 11cm or larger	1/ea
<b>Scissor Testing Material, 5.5m x 14cm Latex</b>		
SH39973	Yellow, Thin - for scissors 10.5cm or smaller	1/ea
SH39974	Red, Medium - for scissors 11cm or larger	1/ea
SH42405	Scissor testing material holder	1/ea



## See every detail

High quality metal arm for optimal flexibility and durability. 60 ultra-bright daylight™ LEDs for shadowfree lighting with 2 brightness levels for clear viewing. Perfect for attaching to packing tables and benches to view instruments as they're being cleaned

- Provides 1.75x magnification
- Large magnification lens (17.5cm diameter)

Order #	Description	UOM
MEK DAYA25130	Daylight Omega 7 magnifying lamp	1/ea

## Hand and nail scrub brush

The hand and nail scrub brushes have a versatile design that's ideal for multipurpose cleaning applications including scrubbing hands

Order #	Description	UOM
SH27000	Hand and nail brush 104mm, Non-Autoclavable	3/pk
SH105608	Hand and nail brush 127mm	3/pk



SH27000



SH105608

## Purekeys flat medical keyboard



Designed for use in all medical environments, combining maximum infection control with excellent tactile key response. Can be used in CSSDs on packing tables and benches

- Completely flat top
- Easy to clean with all hospital grade disinfectants
- Plug and Play compatible with every Windows-based PC

Order #	Description
K-30000100	Medical keyboard compact, wired
K-40000100	Medical keyboard compact fixed angle, wired
K-40100100	Medical keyboard compact with touchpad, fixed angle, wired
K-50500100	Medical Keyboard full size, wired

## Autoclave gloves & mitts

Specially designed autoclave gloves and mitts offer superior protection when handling hot instrument trays and containers.

The mitts with gauntlet-style cuffs provide up-to-the elbow protection. Made of soft terrycloth cotton, the mitts absorb moisture and provide long-lasting comfort for frequent wear.

Completely reusable and can be laundered.

- One-size-fits-all
- Packaged in pairs

Order #	Description	UOM
<b>Single-Ended Cleanable 63.5cm</b>		
SH19586	Autoclave gloves, 28cm gauntlet cuff	1/pair
SH19587	Autoclave mitts, 23cm gauntlet Cuff	1/pair



## Vicolab® silicone gloves



Silicone gloves that are design optimised for the demands of the sterilisation department.

- Thermo - resistant from - 60°C to + 300°C
- 100% autoclavable, sterilisable and washable
- CE approved silicone material

Available in six colours: light blue, dark blue, red, green, black, white  
Max Planck Institute of Germany Quality of certification

Order #	Description	UOM	Order #	Description	UOM
<b>Silicone autoclave gloves</b>			<b>Silicone autoclave gloves</b>		
SHV534671-LB	Light blue	1/ea	SHV534671-GR	Green	1/ea
SHV534671-DB	Dark blue	1/ea	SHV534671-BK	Black	1/ea
SHV534671-RD	Red	1/ea	SHV534671-WH	White	1/ea

# CSSD Personnel Protection



## 3D curved design for “better” shield protection

If your head has a curve so should your face shield, right? Right! The BetterShield Full Face Shield is a contoured face shield that provides upper-, middle-, and lower-face protection, comfortably. Vents let hot air escape and the unique material won't fog up, under any circumstance. The elastic headband provides a secure fit

- Innovative design protects entire face, including sides and bottom
- Distortion-free optics
- Fits easily over glasses, goggles, and N95 masks

Order #	Description	UOM
<b>BetterShield™ Full Face Shield</b>		
SH11586	Small Pack	25/pk
SH17296	Large Pack	100/pk
SH18028	Foam BetterShield™ Full Face Shield	10/pk

## Splash shield disposable

Face Protection System with maximum face protection.

- Anti-fog
- Optically clear
- Lightweight for comfort with new soft strap
- 33cm impenetrable foam barrier for wide coverage

Order #	Description	UOM
<b>Splash Shield</b>		
SFS4505	Full Face - 33cm wide x 15cm long lens	24/pk
SFS4511	Short Face - 33cm wide x 24cm long lens	24/pk



**SFS4505**  
Full Face Splash Shield



**SFS4511**  
Short Face Splash Shield


## Glove box holders




PEG-S-GBH-ABS



PEG-TR-GBH-ABS



PEG-GLOVEH



PEG-CRTA-GLOVBX-3

Wire and solid glove box holders to hold single box or 3 boxes.

- Easy to clean and wash
- Multi mounting options including wall, carts, beds, cabinets
- Quick and easy to switch empty box with new box
- Promotes infection control practices

Solid plastic box holder:

- Fits large variety of boxes (maximum box size: 8x24x14.5cm)
- Comes with S, M, L stickers to affix to front and identify size of gloves

Wire box holder:

- Reduces dust downs

Order #	Description	UOM
<b>Glove box holder solid plastic</b>		
PEG-S-GBH-ABS	Single - 11x26x15cm	1/ea
PEG-TR-GBH-ABS	Triple set - 11x26x46cm	1/ea
<b>Glove box holder wire</b>		
PEG-GLOVEH	Glove holder for boxes 16.4x24x9.2cm	1/ea
PEG-CRTA-GLOVBX-3	Glove holder - 3 tiers for each box 16.4x24x9.2cm	1/ea



# CSSD Personnel Protection



## Eye and face shields

Lightweight and comfortable, the frames are reusable and come in mixed bright colours.

The lenses are disposable - simply pull off a dirty lens and replace it with a new one.

Eye and Face Shields fit over eyeglasses comfortably.

- Distortion-free optical lenses
- Reusable frames come in a variety of colours
- Easily dispose and replace dirty lenses

Snap on.  
Snap off



Lens attachment or removal  
is as easy as 1-2-3!



Order #	Description	UOM	Image
S1025	Clear Choice Eye Shields - FRAMES & LENSES - 10 frames / 25 lenses	Pack 10/25	
SLF100	Clear Choice Eye Shields - FRAMES & LENSES WITH PROTECTIVE FILM - assorted colours - 100 frames / 100 lenses	Pack 100/100	
SD100	Clear Choice Eye Shields - LENSES ONLY - green film - dispenser pack	Pack 100	
SDN100	Clear Choice Eye Shields - LENSES ONLY - clear - dispenser pack	Pack 100	
SDC100	Clear Choice Eye Shields - LENSES ONLY - double green film - dispenser pack	Pack 100	
SH15161	Clear Choice Eye Shield - FRAMES ONLY - black only	Pack 100	
SBF100	Clear Choice Eye Shields - FRAMES ONLY - folding assorted colours	Pack 100	

Order #	Description	UOM	Image
SF1010	Clear Choice Face Shields - FRAMES & LENSES - 10 frames / 10 lenses	Pack 10/20	
SF2510	Clear Choice Face Shields - LENSES ONLY - green protective film	Pack 25	
SFS4510	Clear Choice Face Shields - LENSES ONLY - clear lenses	Pack 100	
SF2511	Clear Choice Face Shield - FRAMES ONLY - black only with adjustable band	Pack 25	
SH17411	Clear Choice Foam Face Shield - FRAME & LENSE - 19cm long lens	Pack 40	
SFF3510	Clear Choice Foam Face Shield - FRAME & LENSE - 24cm long lens	Pack 40	
S9001	Clear Choice Eye Shield & Frame - DISPENSER - wire	Each	

\* green film peels off before use

# NSC for CSSD Departments

National Surgical Corporation (NSC) has worked with various CSSD departments around Australia and understands the requirements of CSSD Managers and Technicians.

The NSC range of CSSD cleaning brushes and equipment is selected from the Sharn USA range and other key suppliers who are leading providers of sterile processing and instrument care supplies. We aim to supply only the highest quality CSSD equipment and cleaning brushes.

The complete NSC range of products and equipment specific for CSSD and sterilisation departments includes:

- Cleaning brushes
- Protection equipment
- Instrument tip protection
- CSSD cleaning accessories
- Steriliser maintenance products
- Sterilisation trays and baskets
- CSSD height adjustable packing tables
- Shelving and storage solutions - high density modular storage, chrome wire shelving

Please contact us to arrange a discussion with an NSC Account Manager to identify your specific CSSD requirements.

## CSSD height adjustable packing tables

*high quality, customised, adjust height as required*



Designed specially for any hospital, clinic or day surgery sterilisation/ CSSD department

- Australian designed and manufactured
- Various sizes available up to 2.3m
- Comes standard with castors and electric linak lifting legs
- Press button to raise or lower table to height (750mm - 1400mm) required by the user, to reduce the risk of staff OH&S issues
- Variety of options including overhead magnifying lamp, under or over bench shelves, roll holders and various others

## Pegasus high density modular storage systems

*essential for any CSSD department to safely store sterile trays*



- High quality interchangeable module system
- Custom designed to suit the storage area
- Can be used in any CSSD store room, control & picking area and general areas
- Improve organisation and cleanliness
- Ability to move modules around for stock rotation
- Install with wire shelves for holding sterile trays (up to 13/bay)
- Baskets to be viewed on an angle without falling
- Perfect system essential for CSSD departments



## Chrome wire static shelving units

*high quality, medical grade chrome wire shelving which supports heavy and bulky loading*



- High quality chrome wire shelving
- Designed specifically for the Australian medical market
- Heavy duty commercial grade wire shelving with support bridges under each shelf
- Heavy and bulky loading of up to 150kg per shelf, 500kg per unit
- Ideal for CSSD - sterile store rooms, control & picking rooms, loan store rooms, deboxing rooms
- Custom designed to suit the storage area

## Chrome wire trolley units

*high quality, heavy duty, chrome wire trolley which enables items to be easily moved around areas*



- High quality chrome wire trolleys with heavy duty castors
- Customised to suit your requirements
- Heavy duty commercial grade wire shelving with support bridges under each shelf
- Heavy and bulky loading of up to 150kg per shelf, 400kg per unit
- Allow safe movement of goods around the hospital
- Ideal for CSSD - sterile store rooms, control & picking rooms, loan store rooms, deboxing rooms

## NSC top track Compactors

*high quality, custom made top track compactors, ideal for storing goods in areas with space restrictions*



- Custom designed and made top track compactors
- Increase storage space by up to 40%
- Ball bearings ensures streamlined movement of mobile bays
- Can be installed in CSSD areas

### **NSC: providing a complete shelving & storage solution**

Our Account Managers will visit your site to discuss your individual requirements and provide a solution to maximise use of your space.

NSC can provide a solution utilising our entire range of shelving and storage solutions and provide detailed CAD designed drawings to demonstrate the correct scale and size for the shelving site.

We have worked in the healthcare industry for over 30 years and understand the requirements and demands of a CSSD department.

# Quick Reference Guide

Unsure what brush you may need? Use the table below to determine what you require based on the diameter

Required	Product Code/ Page (#)
1mm	SH30618 (6)
1.5mm	SH18898 (8), SH21474 (8), SH21479 (8), SH18901 (8)
1.67mm	SH54759 (6)
2mm	SH22172 (10), SH45-542SS (5), SH22173 (10), SH30620 (6), SH54844 (6), SH62387 (9), SH62427 (9), SH62395 (9), SH62436 (9), SH62405 (9), SH62447 (9), SH62415 (9), SH62457 (9), SH105602 (11)
2.3mm/2.4mm	2.3mm: SH30621 (6); 2.4mm: SH18897 (8), SH18900 (8), SH18902 (8), SH21478 (8)
2.5mm/2.7mm	2.5mm: SH45-550SS (5), BR-IM-3 (5), SH34317 (6); 2.7mm: SH30622 (6)
3mm	SH22171 (10), SH22170 (10), SH45-543SS (5), SH54845 (6), SH30623 (6), SH62388 (9), SH62428 (9), SH62396 (9), SH62437 (9), SH62406 (9), SH62448 (9), SH62416 (9), SH62458 (9), SH105603 (11)
3.2mm/3.3mm	3.2mm: SH18896 (8), SH21477 (8), SH21729 (8), SH21730 (8); 3.3mm: SH30624 (6)
3.7mm	SH30625 (6)
4mm	SH45-544-2SS (5), BR-IM-4 (5), SH30626 (6), SH34306 (6), SH62389 (9), SH62430 (9), SH62397 (9), SH62439 (9), SH62407 (9), SH62449 (9), SH62417 (9), SH62459 (9)
4.2mm	SH18899 (8), SH21731 (8)
5mm	BR-45-447 (7), SH45-539SS (5), SH45-545-2SS (5), SH45-551SS (5), BR-45-667 (7), SH20659 (6), SH34310 (6), SH34318 (6), SH18621 (11), SH22172 (10), SH22167 (10), SH22166 (10), SH22162 (10), SH22164 (10), SH22163 (10), SH22165 (10), SH27029 (11), SH27032 (11), SH27033 (11), BR-IM-5 (5), SH62390 (9), SH62431 (9), SH62398 (9), SH62440 (9), SH62408 (9), SH62450 (9), SH62418 (9), SH62463 (9), SH62419 (9), SH62464 (9)
6mm	SH45-546-2SS (5), BR-IM-6 (5), SH34308 (6), SH62391 (9), SH62432 (9), SH62399 (9), SH62441 (9), SH62409 (9), SH62451 (9), SH62420 (9), SH62465 (9)
7mm	SH45-540SS (5), SH45-552SS (5), BR-45-665 (7), BR-IM-7 (5), SH27029 (11), SH34312 (6), SH34319 (6)
8mm	SH62392 (9), SH62433 (9), SH62400 (9), SH62442 (9), SH62410 (9), SH62452 (9), SH62421 (9), SH62466 (9) SH18443 (11)
10mm	BR-45-448 (7), SH45-541SS (5), SH45-553SS (5), BR-45-668 (7), SH34309 (6), SH34313 (6), SH34320 (6), SH62393 (9), SH62434 (9), SH62401 (9), SH62443 (9), SH62412 (9), SH62453 (9), SH62422 (9), SH62467 (9)
12mm	SH27032 (11), SH27033 (11), SH45-541-4SS (5), SH34314 (6), SH34321 (6), SH45-541-4SS (5)
12.5mm	SH62394 (9), SH62435 (9), SH62402 (9), SH62444 (9), SH62412 (9), SH62454 (9), SH62424 (9), SH62468 (9)
15mm	BR-45-449 (7), BR-45-669 (7), SH62403 (9), SH62445 (9), SH62413 (9), SH62455 (9), SH62425 (9), SH62469 (9)
20mm	SH62404 (9), SH62446 (9), SH62414 (9), SH62456 (9), SH62426 (9), SH62470 (9)
Sets/Kits	SH103630 (6), SH103631 (6), SH105602 (11), SH105603 (11)
Double ended	SH45-755SS (7), SH22163 (10), SH22165 (10), SH22169 (10), SH27029 (11), SH27032 (11), SH27033 (11), SH54848 (4), SH54849 (4)



# Index

Code	PG #
6948	13
76.16.03	21
76.16.04	21
76.16.05	21
ACC-5010	18
ACC-5011	18
ACC-5012	18
ACC-5013	18
ACC-5014	18
ACC-5015	22
ACC-5015-SP	22
ACC-5016	22
BR-45-447	7
BR-45-448	7
BR-45-449	7
BR-45-667	7
BR-45-668	7
BR-45-669	7
BR-IM-10	5
BR-IM-3	5
BR-IM-4	5
BR-IM-5	5
BR-IM-6	5
BR-IM-7	5
CP-01	23
CSSD-SIR-1	21
CSSD-SIR-2	21
K-30000100	27
K-40000100	27
K-40100100	27
K-50500100	27
L4060-HR	15
LT-100	19
LT-110	19
LT-112	19
LT-116	19
LT-200	19
LT-210	19
LT-300	19
LT-310	19
LT-400	19
LT-450	19
LT-500	19
LT-550	19
LT-600	19
LT-700	19
LT-720	19
LT-800	19
LT-900	19
LT-930	19
MEK DAYA25130	26
PEG-CRTA-GLOVBX-3	28
PEG-GLOVEH	28
PEG-S-GBH-ABS	28

Code	PG #
PEG-TR-GBH-ABS	28
PP-JPC-1086	13
PP-JPC-1086-B50	13
PP-JPC-3NS	13
PP-JPC-3NSB	13
Rack P	5
S1025	29
S9001	29
SBF100	29
SD100	29
SDC100	29
SDN100	29
SF1010	29
SF2510	29
SF2511	29
SFF3510	29
SFS4505	28
SFS4510	29
SFS4511	28
SH103630	6
SH103631	6
SH103697	7
SH104764	25
SH105500	25
SH105503	25
SH105602	11
SH105603	11
SH105604	4
SH105605	4
SH105606	4
SH105607	4
SH105608	26
SH111401	18
SH111402	18
SH111403	18
SH111404	18
SH111405	18
SH111406	18
SH111408	18
SH111409	18
SH111410	18
SH111411	18
SH111412	18
SH111413	18
SH111414	18
SH111415	18
SH111416	18
SH111417	18
SH111418	18
SH111419	18
SH111420	18
SH111421	18
SH111422	18
SH111423	18

Code	PG #
SH111424	18
SH111425	18
SH111426	18
SH111427	18
SH111428	18
SH111429	18
SH111430	18
SH111431	18
SH111432	18
SH111433	18
SH111434	18
SH111435	18
SH111437	18
SH111439	18
SH111440	18
SH111441	18
SH111442	18
SH111443	18
SH111444	18
SH111445	18
SH111446	18
SH111447	18
SH111448	18
SH111449	18
SH111450	18
SH11586	28
SH15161	29
SH16603	25
SH16975	23
SH17296	28
SH17411	29
SH18028	28
SH18441	8
SH18442	8
SH18443	11
SH18444	22
SH18445	22
SH18621	11
SH18896	8
SH18897	8
SH18898	8
SH18899	8
SH18900	8
SH18901	8
SH18902	8
SH19586	27
SH19587	27
SH20056	12
SH20195	16
SH20196	16
SH20197	16
SH20198	16
SH20199	16
SH20200	16

Code	PG #
SH20201	16
SH20202	16
SH20203	16
SH20204	16
SH20236	17
SH20237	17
SH20238	17
SH20239	17
SH20240	17
SH20241	17
SH20242	17
SH20243	17
SH20271	15
SH20297	16
SH20298	16
SH20299	16
SH20300	16
SH20301	16
SH20302	17
SH20303	17
SH20304	17
SH20305	17
SH20306	17
SH20307	17
SH20659	6
SH21474	8
SH21477	8
SH21478	8
SH21479	8
SH21650	4
SH21651	4
SH21729	8
SH21730	8
SH21731	8
SH22162	10
SH22163	10
SH22164	10
SH22165	10
SH22166	10
SH22167	10
SH22168	10
SH22169	10
SH22170	10
SH22171	10
SH22172	10
SH22173	10
SH22175	7
SH23489	13
SH23528	13
SH27000	26
SH27029	11
SH27032	11
SH27033	11
SH30618	6

Code	PG #
SH30620	6
SH30621	6
SH30622	6
SH30623	6
SH30624	6
SH30625	6
SH30626	6
SH34084	13
SH34085	13
SH34086	12
SH34087	12
SH34088	12
SH34089	12
SH34306	6
SH34308	6
SH34309	6
SH34310	6
SH34311	6
SH34312	6
SH34313	6
SH34314	6
SH34315	6
SH34316	6
SH34317	6
SH34318	6
SH34319	6
SH34320	6
SH34321	6
SH39973	26
SH39974	26
SH39975	26
SH39976	26
SH41235	14
SH41570	15
SH41571	15
SH42405	26
SH42406-BK	19
SH42406-BL	19
SH42406-GN	19
SH42406-OR	19
SH42406-PK	19
SH42406-RD	19
SH42406-WH	19
SH42407-BK	19
SH42407-BL	19
SH42407-GN	19
SH42407-OR	19
SH42407-PK	19
SH42407-RD	19
SH42407-WH	19
SH42408	19
SH42409	19
SH42410	19
SH42411	19

# Index

Code	PG #
SH42412	19
SH42413	19
SH43609	14
SH43884	14
SH43885	14
SH45-539SS	5
SH45-540SS	5
SH45-541-4SS	5
SH45-541SS	5
SH45-542SS	5
SH45-543SS	5
SH45-544-2SS	5
SH45-545-2SS	5
SH45-546-2SS	5
SH45-550SS	5
SH45-551SS	5
SH45-552SS	5
SH45-553SS	5
SH45-755SS	7
SH54675	23
SH54676	23
SH54677	23
SH54682	16
SH54683	16
SH54684	16
SH54685	16
SH54686	16
SH54687	16
SH54688	16
SH54689	16
SH54690	16
SH54692	16
SH54693	16
SH54694	16
SH54695	16
SH54696	16
SH54697	16
SH54698	16
SH54699	16
SH54702	16
SH54705	17
SH54706	17
SH54707	17
SH54708	17
SH54709	17
SH54712	17
SH54713	17
SH54714	17
SH54716	17
SH54717	17
SH54718	17
SH54720	17
SH54721	17
SH54724	17

Code	PG #
SH54725	17
SH54726	17
SH54727	17
SH54729	17
SH54731	17
SH54732	17
SH54733	17
SH54734	17
SH54735	17
SH54736	17
SH54737	17
SH54738	17
SH54739	17
SH54740	17
SH54741	17
SH54742	17
SH54743	17
SH54744	17
SH54745	17
SH54746	17
SH54747	17
SH54748	17
SH54749	16
SH54753	16
SH54754	16
SH54755	16
SH54756	16
SH54757	16
SH54758	16
SH54759	6
SH54844	6
SH54845	6
SH54848	4
SH54849	4
SH57587-RD	14
SH61826	24
SH62387	9
SH62388	9
SH62389	9
SH62390	9
SH62391	9
SH62392	9
SH62393	9
SH62394	9
SH62395	9
SH62396	9
SH62397	9
SH62398	9
SH62399	9
SH62400	9
SH62401	9
SH62402	9
SH62403	9
SH62404	9

Code	PG #
SH62405	9
SH62406	9
SH62407	9
SH62408	9
SH62409	9
SH62410	9
SH62411	9
SH62412	9
SH62413	9
SH62414	9
SH62415	9
SH62416	9
SH62417	9
SH62418	9
SH62419	9
SH62420	9
SH62421	9
SH62422	9
SH62424	9
SH62425	9
SH62426	9
SH62427	9
SH62428	9
SH62430	9
SH62431	9
SH62432	9
SH62433	9
SH62434	9
SH62435	9
SH62436	9
SH62437	9
SH62439	9
SH62440	9
SH62441	9
SH62442	9
SH62443	9
SH62444	9
SH62445	9
SH62446	9
SH62447	9
SH62448	9
SH62449	9
SH62450	9
SH62451	9
SH62452	9
SH62453	9
SH62454	9
SH62455	9
SH62456	9
SH62457	9
SH62458	9
SH62459	9
SH62463	9
SH62464	9

Code	PG #
SH62465	9
SH62466	9
SH62467	9
SH62468	9
SH62469	9
SH62470	9
SH79494	18
SH87340	14
SH90861	24
SHV534671-BK	27
SHV534671-DB	27
SHV534671-GR	27
SHV534671-LB	27
SHV534671-RD	27
SHV534671-WH	27
SHV558908	22
SHV569104	22
SHV578623	23
SHV588690	22
SIMAT-250X255	22
SIMAT-250X420	22
SLF100	29
SW1-2000	21
SW2-2000	21
SWO-2000	21
T104060-AM	15
T104060-ND-HR	15
T204060-AM	15

**Contact your NSC Account Manager for further information  
and pricing of your CSSD requirements**

**1800 138 138 • [info@nationalsurgical.com.au](mailto:info@nationalsurgical.com.au)**



# Let NSC Support You

**BE SUPPORTED BY AN  
AUSTRALIAN COMPANY**

*PROVIDING VALUE BEYOND THE PRODUCT*

## What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

## Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.

## Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has an extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!

## Our range includes:




Wound Drains




Electrosurgical Range



OT & General Consumables



Surgical Instruments




Laparoscopic Instruments



Medichill Ice & Heat Packs



Protection Products



Cleaning & Sterilisation



Shelving & Storage



Pegasus High Density Storage

## CUSTOMER SERVICE

Contact your National Surgical Corporation Account Manager for more information

**Phone: 1800 138 138 (free call) Fax: 1800 138 140 (free fax)**

**[info@nationalsurgical.com.au](mailto:info@nationalsurgical.com.au)**

**[www.nationalsurgical.com.au](http://www.nationalsurgical.com.au)**