

NATIONAL SURGICAL CORPORATION

Silicone instrument mats



National Surgical Corporation
is certified to
AS/NZS ISO 9001:2015

A proudly founded and owned Australian company - providing value to our customers beyond the product

Silicone instrument mats

National Surgical Corporation have a range of silicone mats which can be used to protect instruments and trays during the cleaning and sterilisation process

HOLDS INSTRUMENTS FOR STERILISATION

Place delicate and standard instruments on mats to ensure no movement during the sterilisation and cleaning process

DELICATE OR VERY DELICATE INSTRUMENTS

Different sizes available to hold delicate instruments with 3mm, 4mm or 8mm distance between knobs

DIFFERENT SIZES AVAILABLE

Large range of styles and sizes available

FLAT MATS

Flat mats can be used to cover instruments or secure the tray to the base shelf to stop movement of the tray or instruments

CAN FIT IN ANY STERILISATION BASKETS

Ideal for inserting into any type of sterilisation tray, basket or container system. Mats can be cut to size to suit any sized tray

PROTECTS INSTRUMENTS

Silicone mats hold instruments securely to ensure they won't move, fall or accidentally be lost during cleaning and sterilisation

WORK STATION MAT

Large work surface table mat protects and preserves instruments when working on stainless steel benchtops

STERILISATION BASKETS AVAILABLE

Large range of sterilisation baskets, trays, micro mesh baskets and scope baskets also available













Image	Code	Size
	SIMAT-250X255	Silicone mesh non slip mat - 255x250mm
	SIMAT-250X420	Silicone mesh non slip mat - 420x250mm
	SHV588690	Silicone cover mat - 540 x 320mm
	SHV578623	Silicone work surface mat - 2400x1200x3mm



Image	Code	Size	Distance bet. knobs	Knob height	Size of hole	Distance bet. holes
	SH18444	330x460mm - EACH	8mm	10mm	ø 8mm	56mm/40mm
	SH18445	330x460mm - PACK 5	8mm	10mm	ø 8mm	56mm/40mm
	ACC-5015	515x240mm	4mm	15mm	ø 8mm	35mm
	ACC-5016	545x565mm	4mm	19mm	ø 8mm	35mm
	ACC-5015-SP	230x240mm	4mm	15mm	ø 8mm	35mm
	SHV558908	500 x 320mm	3mm	17mm	10mm squares 3mm & ø 11mm circles	10/11mm
	SHV569104	Micro - 500 x 320mm	3mm	17mm	ø 3mm & 8mm	10/35mm



SH558908 Silicone mat



SH569104 Micro mat



Let NSC Support You

**BE SUPPORTED BY AN
AUSTRALIAN COMPANY**

PROVIDING VALUE BEYOND THE PRODUCT

What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.

Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has an extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!

Our range includes:




Wound Drains




Electrosurgical Range



OT & General Consumables



Surgical Instruments



Laparoscopic Instruments




Medichill Ice & Heat Packs



Protection Products



Cleaning & Sterilisation



Shelving & Storage



Pegasus High Density Storage

CUSTOMER SERVICE

Contact your National Surgical Corporation Account Manager for more information

Phone: 1800 138 138 (free call) Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au