

ECHIDNA RACK SYSTEM



NATIONAL SURGICAL CORPORATION

Australian distributor of the Pegasus high density storage range including Echidna racks



**A COMPLETELY SEALED AND CLEANABLE ALUMINIUM RACKING SYSTEM
IDEAL FOR YOUR OPERATING THEATRE AND CENTRAL STERILISATION UNIT**


UNIQUE FEATURES

THE ECHIDNA RACKING SYSTEM WAS DESIGNED TO PROVIDE THE OPTIMUM SHELVING SOLUTION FOR AREAS THAT REQUIRE THE HIGHEST CLEANING STANDARDS

It is ideal for your operating theatre and central sterilisation environment for storage of sterile trays and packs. Echidna racking can also be used in any other storage area throughout the hospital



SEALED STRUCTURE




The Echidna rack system is a completely sealed system, there are no open holes within the frames and therefore all surfaces can easily be cleaned

LABELLING SYSTEM




Labels can be clicked into the frame and remain in position, even if the tray or basket is removed. Labels can incorporate covers to protect any printed labels



MATERIALS


The entire system is made of powder-coated aluminium; therefore, it can be cleaned with a variety of liquids



100% MODULAR

This completely modular system provides tool free, easy, reconfiguration of basket and shelf levels, both for standard side rails and heavyweight telescopic slides.

Simply click off and on



INSTANT ASSEMBLY

Easily install and modify as required to suit different storage parameters

TELESCOPIC SLIDES


Echidna racks can be equipped with either telescopic slides (for loads over 10kg) or standard plastic side rails with horizontal or angled position

ECHIDNA RACK OPTIONS

THE ECHIDNA RACK IS OFFERED IN TWO ORIENTATIONS, U TYPE AND E TYPE.


The two orientations enable NSC to design a solution to maximum storage capacity within a given area. A U type rack will provide more storage on a given wall while an E type will improve visibility of products and provide a solution for narrow corridors

E TYPE



E type rack single bay:
437 x 690 x 2015mm
Each additional bay:
437 x 637 x 2004mm

U TYPE



U type rack single bay:
637 x 490 x 2015mm
Each additional bay:
637 x 446 x 2004mm

MODULES


The Echidna racks can be installed with various modules, baskets and trays. Use dividers and label holders to further enhance organisation of stored items




**Stainless steel
wire basket 10cm**




**Stainless steel
wire basket 20cm**




**Flat shelf
60x40cm**




**Wire instrument shelf
Chrome or Stainless Steel***



**Chrome wire
basket 10cm**



**Chrome wire
basket 20cm**



**ABS basket
cream 10cm**



**ABS basket
cream 20cm**



**Antimicrobial
green tray 10cm**



**Antimicrobial
green tray 20cm**




**High resistance
yellow tray 10cm**



**High resistance
yellow tray 20cm**




**ABS tray
cream 5cm**




**Non slip silicone mats
42x25cm, 25x25.5cm**



**Stainless steel
U flat container shelf***



**Dividers and
label holders**



** only available on U type racks*



ECHIDNA CARTS

The Echidna cart is lighter than most steel carts, flexible and easily manoeuvrable. It is built with an aluminium structure, HPL top working surface and dual push handles

THE CART WILL MAINTAIN ITS PERFECT AESTHETIC LOOK OVER THE YEARS

Hospital theatre • CSSD • Set up trolley • Store rooms • Emergency Dept • Wards • Day Surgeries • Clinics

- ✓ **100% Aluminium structure**
- ✓ **Lightweight & versatile**
- ✓ **Single or double carts**
- ✓ **Sealed structure – no growth of bacteria inside the profiles**
- ✓ **99cm working surface**
- ✓ **Easy to clean**



Product Code	Description	Dimensions (DxWxH)
HYG-CRT98-E	Echidna cart single E	740 x 490 x 1060mm
HYG-CRT98-UU	Echidna cart double U	690 x 990 x 1060mm

The Echidna open cart is ideal for handling and transporting sterile trays, packs, consumables and containers between CSSD and theatre

This lightweight cart can be used as a clean or dirty cart combining a range of modules including flat wire shelves, high resistance trays, or ABS cream baskets. Available in single or double carts



Echidna open storage carts with yellow high resistance modules or antimicrobial green trays can be used for transportation of used, contaminated instruments, trays and equipment from theatre back to CSSD



Echidna open storage carts with sterile trays and packs can be used in sterile environments and “clean” corridors between CSSD, CSSD store rooms or theatre set up

Product Code	Description	Dimensions (DxWxH)
HYG-CRT160-1-OP	Echidna open cart single	690 x 620 x 1600mm
HYG-CRT160-2-OP	Echidna open cart double	690 x 1063 x 1600mm

** width includes handles*



Let NSC Support You

**BE SUPPORTED BY AN
AUSTRALIAN COMPANY**

PROVIDING VALUE BEYOND THE PRODUCT

What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.


Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has an extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!


Our range includes:




Wound Drains




Electrosurgical Range



OT & General Consumables



Surgical Instruments



Laparoscopic Instruments




Medichill Ice & Heat Packs



Protection Products



Cleaning & Sterilisation



Shelving & Storage



Pegasus High Density Storage

CUSTOMER SERVICE

Contact your National Surgical Corporation Account Manager for more information

Phone: 1800 138 138 (free call) Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au