




GREENMED

Preserving value, redesigning waste.



**100%
RECYCLED &
RECYCLABLE**



GREENGUARD PROTECTIVE CORNER

Product code: GM00300 (box of 300)

- ♻️ Developed in consultation with hospital sterilisation departments
- ♻️ Elevates the tray to allow optimised circulation of sterilant
- ♻️ Manufactured from 100% recycled #5 Polypropylene
- ♻️ Fully recyclable via co-mingled recycling or the GreenMed *Keep Me In The Loop*® Program
- ♻️ Single use guarantees quality every time

Proudly distributed by



**NATIONAL
SURGICAL
CORPORATION**

- ☎ 1800 138 138
- ✉ info@nationalsurgical.com.au
- 🌐 www.nationalsurgical.com.au