



Electrosurgical Plates

National Surgical Corporation Electrosurgical Plates: cost effective, high quality, disposable

HIGH QUALITY ADVANCED GEL ADHESIVE

Provides consistent thermal distribution, conforms to the patient, leaves no residue and is easily removed without damaging friable skin

COMPATIBILITY

All NSC electrosurgical plates/return electrodes are split plates, compatible with all modern electrosurgical generators and monitoring systems

ECONOMICAL CABLE OPTIONS

Available with either an integrated disposable cable or separate reusable cable

SINGLE USE

Disposable to allow for quick and efficient clean up after surgery

VARIOUS SIZES AVAILABLE

Available in two adult styles, with or without cable

ELECTROSURGICAL ACCESSORIES

Electrosurgical pencils, electrode tips, scratch pads and smoke evacuator pencils are also available



National Surgical Corporation
is certified to
AS/NZS ISO 9001:2015



Electrosurgical Plates



Image	Product Code	Product Description	Pad Size	Unit of Purchase
ADULT				
	ESE-SY-MW2	Return electrode pad, single use, multi element with cable	202mm x 107mm	Box 40
	ESE-SY-M2	Return electrode pad, single use, multi element without cable	202mm x 107mm	Box 100
	ESE-SY-MW4	Return electrode pad, single use, multi element with cable	143mm x 185mm	Box 40
	ESE-SY-M4	Return electrode pad, single use, multi element without cable	143mm x 185mm	Box 100
CABLES - REUSABLE				
	ESE-SY-C1	Return electrode cable 3m standard plug	For use with pads without packaged cables	EACH
	ESE-SY-C1-5	Return electrode cable, 5m standard plug		EACH

CUSTOMER SERVICE
 Phone: 1800 138 138 (free call)
 Fax: 1800 138 140 (free fax)
 info@nationalsurgical.com.au
 www.nationalsurgical.com.au

Contact your National Surgical Corporation Account Manager for more information



NATIONAL SURGICAL CORPORATION

A proudly founded and owned Australian company - providing value to our customers beyond the product