

# Bipolar forceps

# Bipolar forceps

Bipolar forceps are specialised surgical instruments used in surgical procedures, especially in delicate or precise surgeries. They work by applying electrical current through the tips of the forceps to coagulate tissue, stop bleeding, or cut through tissue. The main purpose of bipolar forceps is to control bleeding by cauterising blood vessels during surgery.

Bipolar forceps are used in various surgical procedures including neurosurgery, laparoscopic surgery, ENT (Ear, Nose, and Throat) surgery, general surgery and gynecological surgery.

NSC bipolar forceps are sterile, single use and provided in Box 25

## AMERICAN BIPOLAR FORCEPS WITH CABLE

### McPherson straight

9cm long, 3m cable



Code	Tip size
B-FR-700-0902C	0.25mm
B-FR-700-0905C	0.5mm
B-FR-700-0910C	1.0mm
B-FR-700-0915C	1.5mm
B-FR-700-0920C	2.0mm

### Iris straight

11cm long, 3m cable



Code	Tip size
B-FR-710-1102C	0.25mm
B-FR-710-1105C	0.5mm
B-FR-710-1110C	1.0mm
B-FR-710-1115C	1.5mm
B-FR-710-1120C	2.0mm

### Adson

12cm long, 3m cable



Code	Tip size
B-FR-720-1202C	0.25mm
B-FR-720-1205C	0.5mm
B-FR-720-1210C	1.0mm
B-FR-720-1215C	1.5mm
B-FR-720-1220C	2.0mm

## Cushing

3m cable



Code - 15.2cm	Tip size	Code - 18cm	Tip size	Code - 20cm	Tip size
B-FR-730-1502C	0.25mm	B-FR-730-1802C	0.25mm	B-FR-730-2002C	0.25mm
B-FR-730-1505C	0.5mm	B-FR-730-1805C	0.5mm	B-FR-730-2005C	0.5mm
B-FR-730-1510C	1.0mm	B-FR-730-1810C	1.0mm	B-FR-730-2010C	1.0mm
B-FR-730-1515C	1.5mm	B-FR-730-1815C	1.5mm	B-FR-730-2015C	1.5mm
B-FR-730-1520C	2.0mm	B-FR-730-1820C	2.0mm	B-FR-730-2020C	2.0mm


## Cushing angled

3m cable



Code - 15.2cm	Tip size	Code - 18cm	Tip size	Code - 20cm	Tip size
B-FR-731-1502C	0.25mm	B-FR-731-1802C	0.25mm	B-FR-731-2002C	0.25mm
B-FR-731-1505C	0.5mm	B-FR-731-1805C	0.5mm	B-FR-731-2005C	0.5mm
B-FR-731-1510C	1.0mm	B-FR-731-1810C	1.0mm	B-FR-731-2010C	1.0mm
B-FR-731-1515C	1.5mm	B-FR-731-1815C	1.5mm	B-FR-731-2015C	1.5mm
B-FR-731-1520C	2.0mm	B-FR-731-1820C	2.0mm	B-FR-731-2020C	2.0mm

### Tip sizes



Bipolar forceps feature a 2 pin connector

# Bipolar forceps

## AMERICAN BIPOLAR FORCEPS WITH CABLE

### Jeweler #4

11cm long, 3m cable



Code	Tip size
B-FR-740-1102C	0.25mm
B-FR-740-1105C	0.5mm
B-FR-740-1110C	1.0mm
B-FR-740-1115C	1.5mm
B-FR-740-1120C	2.0mm

### Jeweler #5

11.4cm long, 3m cable



Code	Tip size
B-FR-750-1102C	0.25mm
B-FR-750-1105C	0.5mm
B-FR-750-1110C	1.0mm
B-FR-750-1115C	1.5mm
B-FR-750-1120C	2.0mm

### Gerald

18cm long, 3m cable



Code	Tip size
B-FR-780-1802C	0.25mm
B-FR-780-1805C	0.5mm
B-FR-780-1810C	1.0mm
B-FR-780-1815C	1.5mm
B-FR-780-1820C	2.0mm

### Hardy bayonet

20cm long, 3m cable



Code	Tip size
B-FR-800-2002C	0.25mm
B-FR-800-2005C	0.5mm
B-FR-800-2010C	1.0mm
B-FR-800-2015C	1.5mm
B-FR-800-2020C	2.0mm

## Jansen bayonet

3m cable



Code - 15.2cm	Tip size	Code - 18cm	Tip size	Code - 20cm	Tip size
B-FR-760-1502C	0.25mm	B-FR-760-1802C	0.25mm	B-FR-760-2002C	0.25mm
B-FR-760-1505C	0.5mm	B-FR-760-1805C	0.5mm	B-FR-760-2005C	0.5mm
B-FR-760-1510C	1.0mm	B-FR-760-1810C	1.0mm	B-FR-760-2010C	1.0mm
B-FR-760-1515C	1.5mm	B-FR-760-1815C	1.5mm	B-FR-760-2015C	1.5mm
B-FR-760-1520C	2.0mm	B-FR-760-1820C	2.0mm	B-FR-760-2020C	2.0mm

## Semkin

14cm long, 3m cable



Code	Tip size
B-FR-810-1402C	0.25mm
B-FR-810-1405C	0.5mm
B-FR-810-1410C	1.0mm
B-FR-810-1415C	1.5mm
B-FR-810-1420C	2.0mm


## Bipolar cable

For use with REUSABLE bipolar forceps



Code	Dimensions
TAY-L0020-B	3 metres, 2 pin, reusable
AL-SW46211	3 metres, 2 pin, sterile, single use
AL-SW46200A	3 metres, 2 pin, reusable

### Tip sizes



Bipolar forceps feature a 2 pin connector



# Bipolar forceps

## EUROPEAN BIPOLAR FORCEPS WITH CABLE

### Jeweler #4 straight

10.8cm long, 3m cable



Code	Tip size
B-FR-510-1102C	0.25mm
B-FR-510-1105C	0.5mm
B-FR-510-1110C	1.0mm
B-FR-510-1115C	1.5mm
B-FR-510-1120C	2.0mm

### Jeweler #4 curved

10.8cm long, 3m cable



Code	Tip size
B-FR-520-1102C	0.25mm
B-FR-520-1105C	0.5mm
B-FR-520-1110C	1.0mm
B-FR-520-1115C	1.5mm
B-FR-520-1120C	2.0mm

### Jeweler #7 curved

10.8cm long, 3m cable



Code	Tip size
B-FR-531-1102C	0.25mm
B-FR-531-1105C	0.5mm
B-FR-531-1110C	1.0mm
B-FR-531-1115C	1.5mm
B-FR-531-1120C	2.0mm

## Adson

10.8cm long, 3m cable



Code	Tip size
B-FR-540-1102C	0.25mm
B-FR-540-1105C	0.5mm
B-FR-540-1110C	1.0mm
B-FR-540-1115C	1.5mm
B-FR-540-1120C	2.0mm

## Straight

10.8cm long, 3m cable



Code	Tip size
B-FR-550-1102C	0.25mm
B-FR-550-1105C	0.5mm
B-FR-550-1110C	1.0mm
B-FR-550-1115C	1.5mm
B-FR-550-1120C	2.0mm

## Straight angled

10.8cm long, 3m cable



Code	Tip size
B-FR-551-1102C	0.25mm
B-FR-551-1105C	0.5mm
B-FR-551-1110C	1.0mm
B-FR-551-1115C	1.5mm
B-FR-551-1120C	2.0mm

## Bayonet

10.8cm long, 3m cable



Code	Tip size
B-FR-560-1102C	0.25mm
B-FR-560-1105C	0.5mm
B-FR-560-1110C	1.0mm
B-FR-560-1115C	1.5mm
B-FR-560-1120C	2.0mm



# Let NSC Support You

**BE SUPPORTED BY AN  
AUSTRALIAN COMPANY**

*PROVIDING VALUE BEYOND THE PRODUCT*

## What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

## Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.

## Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has an extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!

## Our range includes:




Wound Drains



Electrosurgical Range



Ot & General Consumables



Surgical Instruments



Patient positioning



Medichill Ice & Heat Packs



Protection Products



Cleaning & Sterilisation



Shelving & Storage



Pegasus High Density Storage

## CUSTOMER SERVICE

Contact your National Surgical Corporation Account Manager for more information

**Phone: 1800 138 138 (free call)**

**[info@nationalsurgical.com.au](mailto:info@nationalsurgical.com.au)**

**[www.nationalsurgical.com.au](http://www.nationalsurgical.com.au)**