



**NATIONAL
SURGICAL
CORPORATION**



*A proudly founded and owned Australian company -
providing value to our customers beyond the product*



POLYCARTS

Medical and crash carts for Healthcare



Polycarts

ENVIRONMENTALLY FRIENDLY

Made from 100% polymer which can be recycled

ANTIMICROBIAL

100% antimicrobial is ideal for hospital environments

WASHABLE

Polymers can be easily cleaned and washed without the risk of rust and corrosion

PROTECTION

360° wrap around bumper protects cart and walls

LIGHTWEIGHT & RIGID

Great rigidity even with its lightweight design

PUSH HANDLE



Easily manoeuvre cart with integrated handle

COLOUR OPTIONS

Cart sides and back available in solid colours (red, blue) or handle colours only

MODULAR

Can be built with various drawer combinations in 5cm, 10cm and 20cm heights



POLYCART OPTIONS

All sizes are DxWxH



PLY36-4449-GG
580x840x920mm



PLY41-44449-GG
580x840x1040mm



PLY46-444449-GG
580x840x1170mm



PLY36-4449-BB
580x840x920mm



PLY41-44449-BB
580x840x1040mm



PLY46-444449-BB
580x840x1170mm



PLY36-4449-RR
580x840x920mm



PLY41-44449-RR
580x840x1040mm



PLY46-444449-RR
580x840x1170mm

Handle colours

All Polycarts can be installed with coloured handles



Accessories for carts

POLYCART SPECIALTY CARTS



Broselow paediatric emergency cart

PLY41-33333344-RR

580x840x1040mm

Features 6 x 5cm drawers, 2 x 10cm drawers

Each drawer is designed to hold the emergency equipment necessary to perform resuscitation on a child in a weight range designated by colour. The colour coded draws correspond to the colour zones of the Broselow tape to allow for fast and easy identification during an emergency.

Incorporate with individual drawer locks or full external bar.



Polycart medication cart

PLY41-33C39-BB

580x840x1040mm

Features 1 x 5cm drawers, 3 levels of cassettes with holder, 1 x 20cm drawers

Individual medication cassettes can be stocked for individual patients and removed as needed. Complete medication cassette holder can be removed with all cassettes for integration with other systems or to restock cassettes.

Features electronic lock for additional security of medications.



Polycart IT & multi-purpose cart

PLY41-44449-BB-PO

580x840x1040mm

Features 4 x 10cm drawers, 1 x 20cm drawers, pull out side shelf

Retractable side shelf is pushed to open which allows additional workspace for writing notes and patient information. Monitor screen and various electronic accessories allow laptop, Ipads, VESA screens and other devices to be used and displayed at various heights.



Polycart crash/emergency cart

PLY41-44449-RR

580x840x1040mm

Features 4 x 10cm drawers, 1 x 20cm drawers, external lock

Red coloured cart is highly visible when emergency occurs. Keep cart fully stocked and external lock allows for goods to stay secure until needed. Combine with defibrillator shelf, IV pole and oxygen tank holder.

LOCK OPTIONS




PLY-KEYLCK
Key lock - standard



PLY-SLCK
Breakaway/seal lock



PLY-CODELCK
Electronic combination
code lock




PLY-EXLK-SET
External lock set for
single drawer
** Fitted at factory*

The breakaway emergency external lock's quick opening and reliability boosts speed in emergency response. With a single pull, all drawers can open at once

Individual locks allow for simultaneous use in multiple emergencies and enables easy and fast restocking.

Fast opening with one pull, and easy to see which drawers were utilised



EXLK-ALU-36/41/46
External lock for
all drawers
** Fitted at factory*




Accessories for carts




PLY-CARBACKB
Cardiac back board
(mounts to rear of cart)
** Mounts fitted at factory*



PLY-SDIFARM
Defibrillator shelf with
swivelling arm



PLY-SBSKT
Multi storage wire basket



CRTA-CATH-HLD
Side catheter box with
flexiglass front panel
Lid - CRTA-CATH-LID



S-GBH-ABS
Glove box holder single
(requires mounting hook
DIN-HKS-GBH-SET)




TR-GBS-ABS
Glove box holder triple
(requires mounting hook
DIN-HKS-GBH-SET)



bottle
not included

DLBH-PP
Holder for disinfection bottles
(requires mounting hook
DIN-HKS-GBH-SET)



HAZ-BIB-H
Holder for sharps/
hazard bins



PLY-TOP-PP-MAT
Top protective polymer mat




CRTA-WSBIN
Waste container with lid




PLY-IVPOLE/ PLY-IV-ADON
Polycart IV pole adjustable
(mounts to rear of cart).
IV pole adaptor used when other
accessories installed on rear




OXBH-KIT
Oxygen tank holder
adjustable







TILTBIN-4
Tilt-Bin cassette 4 sections
(requires tilt bin holder)




TILTBIN-5
Tilt-Bin cassette 5 sections
(requires tilt bin holder)



TILTBIN-6
Tilt-Bin cassette 6 sections
(requires tilt bin holder)



PLY-TBH-2
Tilt-Bin holder/posts for
2 rows of tilt bins



PLY-TBH-3
Tilt-Bin holder/posts for
3 rows of tilt bins



PLY-AIO-HLD
All in one arm



PLY-LTH
Laptop arm holder



PLY-VESA
VESA arm holder



PLY-FSHLF-VF
Rear shelf




PLY-SSHLF-R/L
Stainless steel shelf heavy load
push to open
** Fitted at factory*



PLY-KBTR-VF
Keyboard tray, push to open



PLY-SIDE-TB
Side tilt bins and holder
includes 12 bins




PLY-FRAIL
Rear accessory rails (pair)
*Used for mounting items
on rear of cart for easy reach*



PLY-DWR-LBCOV
Label cover (Pack 60)
slides into handle




PLY-A4PRF-LB
Perforated sheet of labels
printable, slide into
label cover



PLY-BRSLW-SET
Broselow/paediatric
emergency coloured label set

ACCESSORIES MAY REQUIRE ACCESSORY RAILS (1 OR 2) OR POLE SET



PLY-SRAIL
Polycart side accessory rail



PLY-POLE-KIT
Pole kit rear (pair)

- 1 accessory rail required
- 2 accessory rails required
- Pole kit required



Let NSC Support You

**BE SUPPORTED BY AN
AUSTRALIAN COMPANY**

PROVIDING VALUE BEYOND THE PRODUCT

What makes NSC different to our competitors to support you?

- A proudly Australian owned and founded company
- Combining 25 years of expertise in the healthcare industry
- Creating value for our customers beyond the product
- Working with our customers to offer specialised solutions
- Offering a large product range focused on cost effectiveness, patient care and productivity
- Low overheads allows us to provide cost savings to our customers compared to the multinationals
- Coordinated, whole company approach to every customer
- Establishing trust for long term customer relationships

Who are we?

National Surgical Corporation (NSC) is an import manufacturer in the healthcare industry, distributing surgical instruments, consumables, equipment and shelving all major hospitals and healthcare facilities throughout Australia. Our extensive product range can be used by all specialist disciplines and departments within a hospital, day surgery or specialist centre.

Let NSC support you and your department

Don't be afraid to change, NSC will support and assist you throughout every step. NSC has extensive training and inservicing material available and can tailor a program to meet your individual needs so that all personnel are aware of how the product operates. Switching to NSC can save you significantly and we can show you how!


Our range includes:




Wound Drains



Electrosurgical Range



OT & General Consumables



Surgical Instruments



Patient Positioning



Medichill Ice & Heat Packs



Protection Products



Cleaning & Sterilisation



Shelving & Storage



Pegasus High Density Storage

CUSTOMER SERVICE

Contact your National Surgical Corporation Account Manager for more information

Phone: 1800 138 138 (free call)

info@nationalsurgical.com.au

www.nationalsurgical.com.au